



southpaw enterprises, inc.®

Store these instructions along with the enclosed maintenance checklist in a safe place. You may also access them on our website.

Instruction Sheet:

4-in-1 Climbing Steps 150130



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location. The Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST:

- 1 - Set of four climbing steps
- 1 - Instruction Sheet
- 1 - Maintenance Checklist

WORKING LOAD: 200 lbs

MATTING PRECAUTIONS: Always use appropriate matting.

MAINTENANCE CHECKLIST:

All of the following wear points should be checked on a monthly basis and recorded on the checklist enclosed.

- Check for chips or splinters. If chips or splinters are found, sand smooth and refinish with lacquer.
- Check for cracks or loose joints. If cracks or loose joints are found, discontinue use until repaired.

CLEANING INSTRUCTIONS:

Check with your infectious disease department before using these procedures.

- Wipe with a damp cloth using disinfectant or mild soap and water.
- Dry before using.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.