


## custom walker adjustments

*With your arm extended down, handle should be at wrist height.*

*With your hand placed on the hand-grip, your elbow should be flexed 20° - 30°.*



### CAUTION:

Be aware of the following:

- Electrical cords in way
- Liquids on floor
- Throw rugs
- Toys or debris
- Furniture
- Doors and door steps
- Gaps in floor

### WARNING:

1. Use wheels on front legs only or on all four legs with appropriate rear wheel brakes.
2. Do not push one side of the walker ahead and then other side without lifting the walker first.
3. Use Drive accessories only. Use of other manufacturers' accessories can result in improper fit and risk fall and/or injury.
4. Do not use on stairs or escalators.

## lifetime limited warranty

Your Drive brand product is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

In the event of a defect covered by this warranty, we will, at our option, repair or replace the device.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.

deluxe  
trigger release  
folding walker

# walker instructions

easy “pull up” trigger release  
for opening and closing




Item # 10224-2 Shown  
For use with all  
Deluxe Folding, Trigger Release, Folding Walkers.

**WEIGHT CAPACITY: 350 LBS.**


## wheelchair to walker transfer

Making certain wheels on the wheelchair are locked lift your body by pushing down on the armrests and carefully transfer one hand at a time to the walker.



## wheelchair to walker transfer (with lower extremity injury to right leg or foot)

Making certain wheels on the wheelchair are locked, lift your body by pushing down on the armrests and carefully transferring first the left hand, then the right hand to the walker. (If injury to left leg or foot transfer the right hand first, followed by left.)



## three-point gait (for walking with walker)

### Step 1

Align the middle of your foot with the back legs of the walker.



### Step 2

Lift the walker up and place it comfortably in front of you, making sure all four legs end up on the floor.



### Step 2

Carefully, walk forward one step at a time.

