**Owner's Booklet
Manual del Usuario**

FreeStyle

Lite

Blood Glucose Monitoring System
Sistema de Control de Glucosa en la Sangre

 **Abbott**

The FreeStyle Lite® Meter and How It Works

System Check Screen

This display always appears when the meter is turned on. You should check that your meter matches the example exactly every time your meter turns on. **Do not use the meter if the display**

check screen does not exactly match the example. If the screen does not match, the meter may show an incorrect result. Please contact Customer Care.

Display Screen

Displays your test results and other important information.

m (Mode) Button

- Moves to different mode settings
- Scrolls backward
- Silences a reminder alarm

Light Button

Controls the display backlight and test strip port light.

Sample Areas (dark-colored half-circles on test strip)

Apply blood or control solution to one sample area only. Insert strip into the meter printed side up (see image to the right).

Data Port

Downloads test results. Requires a compatible data management system. *For more information, please go to our website or contact Customer Care.*

c (Configure) Button

- Marks a control solution test
- Scrolls forward
- Displays clock
- Silences a reminder alarm

Test Strip Port

Insert the **Top** end of a new FreeStyle Lite® test strip here. The meter powers on when you insert the test strip.


Table of Contents

The FreeStyle Lite® Meter and How It Works	i
Important Information About Your Blood Glucose Monitoring System	1
Intended Use	1
Setting Up the Meter	2
Purpose	2
Set Sound Volume	3
Set Time	4
Set Date	5
Using Meter Lights	7
Testing Your Blood Glucose	8
Purpose	8
Prepare to Test	8
Important Test Strip Information	10
Performing a Blood Glucose Test	11
Understanding Blood Glucose Test Results	13
Low and High Blood Glucose Results	14
Error Codes	17
Troubleshooting	20

FreeStyle® Control Solution	22
Purpose	22
Control Solution Information	22
Out of Range Control Solution Results	23
Performing a Control Solution Test	24
Using Your Meter Memory	27
Purpose	27
Viewing Blood Glucose Averages in Memory	27
Viewing Test Results in Memory	28
Setting Reminder Alarms (optional)	29
Purpose	29
Transferring Test Results To a Computer	30
Maintaining the Meter	31
Replacing the Battery	31
Cleaning the Meter	33
Specifications	34
Definition of Symbols	36

Important Information About Your Blood Glucose Monitoring System

Intended Use

Use outside the body only (*in vitro* diagnostic use).

The FreeStyle Lite® Blood Glucose Monitoring System is intended for use in the quantitative measurement of glucose in capillary whole blood from the finger, upper arm and palm and venous whole blood. It is intended for use by healthcare professionals and people with diabetes mellitus at home as an aid in monitoring the effectiveness of a diabetes control program. It is not intended for the diagnosis of or screening for diabetes mellitus, and it is not intended for use on neonates or arterial blood.

IMPORTANT:

The device must not be used on more than one person including other family members due to the risk of spreading infection. All parts of the device and its accessories are considered biohazardous and can potentially transmit infectious diseases, even after performing the cleaning procedure.^{1,2}

Healthcare professionals handling this system must always wear gloves and should follow the infection control policies and procedures approved by their facility.

Use only FreeStyle Lite test strips and FreeStyle® control solution with this meter. Other products can produce inaccurate results.

¹"FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010) <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>.

²"CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010) <http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

- The expected glucose level for a non-diabetic, non-pregnant fasting adult is under 100 mg/dL.³
- Two hours after meals, levels should be less than 140 mg/dL.³
- Consult your healthcare professional to determine the range that is appropriate for you.

Warning:

- Test on your fingers if you are testing for low blood glucose (hypoglycemia) or if you suffer from hypoglycemia unawareness.
- This system contains small parts that may be dangerous if swallowed.

Caution: Please read all the instructions provided in this owner's booklet and practice the testing procedures before using the FreeStyle Lite® system. Blood glucose monitoring should be done with the guidance of a healthcare professional.

Setting Up the Meter

Purpose

The purpose of this section is to set the meter sound volume, time and date.


1. With the meter off, press and **hold** **(m)** for **3** seconds until **SET** appears.

2. Press **(m)** to bypass setting the optional reminder alarms. You may choose to set these later (see Setting Reminder Alarms section).

³American Diabetes Association Clinical Practice Recommendations: Diagnosis and classification of diabetes mellitus. Diabetes Care 2005;28 (Suppl.1):S37-S42

Set Sound Volume

With the sound volume on, your meter will beep when:



- The test strip is full of blood
- The test result appears on the display
- An error occurs

You may also choose to silence your meter.

3. Set sound volume

If you prefer:	Press C until you see:	Then
No sound		<ul style="list-style-type: none"> • Press m to save
A low beep		
A louder beep		

Set Time

Time Format

Note: Your meter can display either the 12-hour (1:23P) or the 24-hour (13:23) time formats. If you prefer the 12-hour format, there is no “A” for AM. If setting a PM time, continue to press **(C)** until you see the “P” displayed for PM.

Hour

4. Set hour

- Press **(C)** to display correct hour
- Press **(m)** to save

Minutes

5. Set minutes

- Press **(C)** to display correct minutes
- Press **(m)** to save



Time Format

6. Set time format

- Press **(m)** to accept time format or
- Press **(C)** to change
- Press **(m)** to save

Set Date


Note: The date format can display either:

The month will flash first



7. Set month

- Press **(C)** to display correct month
- Press **(m)** to save

8. Set day


- Press **(C)** to display correct day
- Press **(m)** to save

Date Format

9. Set date format

- Press **(m)** to accept date format or
- Press **(c)** to change
- Press **(m)** to save

Year

10. Set year

- Press **(m)** to accept year or
- Press **(c)** to change
- Press **(m)** to save

- ## 11. Press and **hold** **(m)** for **3** seconds until meter turns off. Meter set up is complete.

Using Meter Lights

You can light the display screen or test strip port by pressing  on meter when the meter is turned on.

To light your display backlight, press .

The light turns off:

- When you press  again
- After 10 seconds of inactivity

To light the test strip port after you've inserted the test strip, press .

The light turns off:

- When you press  again
- After 60 seconds of inactivity
- When the test strip is filled
- When the test strip is removed from the meter

Testing Your Blood Glucose

Purpose

The purpose of this section is to produce an accurate blood glucose test result.

Prepare to Test

1. Select your test site.

Note: Avoid moles, veins, bones, and tendons. Bruising may occur at the test site. If bruising occurs, consider selecting another site.

There are differences in testing on fingers versus alternative sites.

Warning:

Test on your fingers if you are testing for low blood glucose (hypoglycemia) or if you suffer from hypoglycemia unawareness.

Check for low blood glucose when:

- You feel symptoms such as weakness, sweating, shakiness, headache, or confusion
- You eat a delayed meal after taking insulin
- You are following your healthcare professional's advice

**IMPORTANT:**

2. Wash your hands and the test site with warm water and soap to ensure accurate results.
 - Thoroughly dry your hands and the test site.
 - Do *not* use lotion or cream on the test site.

3. Check test strip expiration date.

Do *not* use expired test strips; they may produce inaccurate results.

IMPORTANT: Use only FreeStyle Lite® test strips. Other test strips can produce inaccurate results. The test strips are for single use only.


Important Test Strip Information

Please read the FreeStyle Lite® test strip package insert for more information.

- Use the test strips within the meter's operating range: 40° – 104° F (4° – 40° C).
- Store test strips in a cool, dry place between 40° – 86° F (4° – 30° C).
- Store test strips away from direct sunlight and heat.
- Store test strips in their original vial only.
- Clean and dry your hands before removing a strip from the vial.
- Close the vial cap tightly immediately after removing a test strip.
- Do *not* transfer test strips to another container.
- Do *not* store test strips outside the vial.
- Do *not* bend, cut, or alter test strips.
- Do *not* apply blood to both sides of the test strip.
- Do *not* press the test strip against the test site.
- Do *not* scrape the blood onto the test strip.
- Do *not* apply blood to the flat side of the test strip.
- Do *not* apply blood to the test strip when the test strip is out of the meter.
- Do *not* put blood or foreign objects into the meter.

Warning: Drying agents in the test strip vial may be harmful if inhaled or swallowed and may irritate skin and eyes.


Performing a Blood Glucose Test

1. Remove test strip from vial.
2. Insert test strip to turn meter on.


Note: The meter turns off after 2 minutes of inactivity. Remove and reinsert the unused test strip to restart the meter.

If the meter does not turn on, go to Troubleshooting section.

3. Confirm System Check Screen.

- This display appears when you turn on the meter. **Do not use the meter if the display check screen does not exactly match the example. Contact Customer Care.**
 - See The FreeStyle Lite® Meter and How It Works section for more information.
 - If the system check screen does not appear, go to Troubleshooting section.

4. Obtain a blood drop.

When the  and  appears, use the lancing device to obtain a blood drop. (See FreeStyle® lancing device insert for instructions.)


The test strip requires only a 0.3 microliter of blood for accurate results.

5. Apply and **hold** blood to 1 sample area of the test strip until:

- You see short lines moving clockwise on the display or
- You hear a beep

This indicates the test strip obtained enough blood and the meter is checking your glucose level. If a short line does not appear after 5 seconds, the sample may be too small. You can add blood to the same side only for up to 60 seconds. If the test does not start after applying blood drop, go to Troubleshooting section.

(Example only)

6. View result

The result appears on the display when the test is complete. The time that the meter takes to display a result depends on your blood glucose level. Higher glucose levels require more time.

Note: If you see an “**Er1, Er2, Er3, or Er4**” on the display, it is an error code. Go to Error Codes section.

7. Remove test strip to turn meter off.

Discard used lancet and test strip properly.

8. Wash your hands and test site thoroughly with warm water and soap after handling the meter, lancing device or test strips.

Understanding Blood Glucose Test Results

The meter displays results in mg/dL. The unit of measurement is preset. You cannot change this setting.

IMPORTANT: The meter displays results from 20 – 500 mg/dL. Low or high blood glucose results can indicate a potentially serious medical condition.

Low and High Blood Glucose Results

Result	What It Means	What To Do
lower than 60 mg/dL	Low (hypoglycemic)	Repeat the test with a new test strip, even if you do <i>not</i> feel that you have low blood glucose. If your blood glucose result is still not consistent with your symptoms, contact your healthcare professional and follow his or her treatment advice.
higher than 240 mg/dL	High (hyperglycemic)	Repeat the test with a new test strip, even if you do <i>not</i> feel that you have high blood glucose. If your blood glucose result is still not consistent with your symptoms, contact your healthcare professional and follow his or her treatment advice.

Display	What It Means	Symptoms	What To Do
 Appears when result is: lower than 20 mg/dL	Severe low blood glucose	Weakness Sweating Shakiness Headache Confusion	Follow your healthcare professional's advice to treat low blood glucose. If your blood glucose result does not match how you feel, perform a control solution test to confirm that your meter and strips are working properly.

IMPORTANT:

- Severe dehydration may cause false low results. If you believe you are suffering from severe dehydration, contact your healthcare professional immediately.
- Contact your healthcare professional if you have symptoms that do *not* match your test result, and you have followed the instructions in this owner's booklet.

Display	What It Means	Symptoms	What To Do
<p data-bbox="134 107 200 161">HI</p> <p data-bbox="86 184 262 313">Appears when result is: higher than 500 mg/dL</p>	<p data-bbox="312 101 532 163">Severe high blood glucose</p>	<p data-bbox="587 101 783 231">Fatigue Thirst Excess Urination Blurry Vision</p>	<p data-bbox="827 101 1266 339">Follow your healthcare professional's advice to treat high blood glucose. If your blood glucose result does not match how you feel, perform a control solution test to confirm that your meter and strips are working properly.</p>

Error Codes

English

Message	What It Means	What To Do
Er 1	<ul style="list-style-type: none">• The sample is too small• There may be a problem with the test strip• Meter error• Very low blood glucose lower than 20 mg/dL	<ol style="list-style-type: none">1. Follow your healthcare professional's advice to treat low blood glucose if you have symptoms such as:<ul style="list-style-type: none">• Weakness• Shakiness• Sweating• Headache• Confusion2. Perform a control solution test using a new test strip. If the test results are within the range printed on the test strip vial, retest using blood and a new test strip.3. If the control solution result is out of range or the error reappears, contact Customer Care.

Message	What It Means	What To Do
	<ul style="list-style-type: none">• There may be a problem with the test strip• Meter error• Very high blood glucose higher than 500 mg/dL• HIGH control solution was applied when the temperature was too cold. (Applies only to control solution labeled as HIGH).	<ol style="list-style-type: none">1. Follow your healthcare professional's advice to treat high blood glucose if you have symptoms such as:<ul style="list-style-type: none">• Fatigue• Thirst• Excess urination• Blurry vision2. Perform a control solution test using a new test strip. If the test results are within the range printed on the test strip vial, retest using blood and a new test strip.3. If the control solution result is out of range or the error reappears, contact Customer Care.

Message	What It Means	What To Do
[Er3]	<ul style="list-style-type: none">• Incorrect test procedure. For example, putting blood on the test strip before inserting it into the meter or applying blood before  and  appear• There may be a problem with the test strip• Meter error	<ol style="list-style-type: none">1. Wait for  and  to appear before applying blood or control solution.2. Perform a control solution test using a new test strip. If the test results are within the range printed on the test strip vial, retest using blood and a new test strip.3. If the control solution result is out of range or the error reappears, contact Customer Care.
[Er4]	<ul style="list-style-type: none">• There may be a problem with the test strip• Meter error	<ol style="list-style-type: none">1. Perform a control solution test using a new test strip. If the test results are within the range printed on the test strip vial, retest using blood and a new test strip.2. If the control solution result is out of range or the error reappears, contact Customer Care.

Troubleshooting

	What It Means	What To Do
1. The meter does not enter test mode after inserting a test strip.	<ul style="list-style-type: none">• Test strip is inserted upside down, or bottom end in, or not fully inserted into the meter	<ol style="list-style-type: none">1. Insert test strip print side up, top end in.
	<ul style="list-style-type: none">• No battery is installed• Battery is installed incorrectly	<ol style="list-style-type: none">1. Install battery with (+) facing up.
	<ul style="list-style-type: none">• Dead battery	<ol style="list-style-type: none">1. Replace battery. Reset date and time, if necessary.2. If meter still does <i>not</i> enter test mode, contact Customer Care.
	<ul style="list-style-type: none">• Defective test strip or meter• Blood or foreign objects put in the test strip port	<ol style="list-style-type: none">1. Contact Customer Care.

	What It Means	What To Do
2. The test does not start after applying the blood drop.	<ul style="list-style-type: none">Blood drop is too small	<ol style="list-style-type: none">Add blood to the same sample area within 60 seconds, or repeat the test using a new test strip and a larger blood drop.
	<ul style="list-style-type: none">Defective test stripSample applied after meter turns off	<ol style="list-style-type: none">Repeat the test using a new test strip. Wait for  and  to appear before applying blood or control solution.If problem persists, contact Customer Care.
	<ul style="list-style-type: none">Defective meter	<ol style="list-style-type: none">Contact Customer Care.

FreeStyle® Control Solution

Purpose

FreeStyle control solution is a red liquid that contains a fixed amount of glucose.

The purpose of doing a control solution test is:

- To practice testing *without* using your own blood
- To confirm that the meter and test strips work together properly
- When you are unsure of your blood glucose test results

Contact Customer Care for information on how to obtain control solution.

Control Solution Information

- Use only FreeStyle control solution with the meter.
- Replace the cap securely on the bottle immediately after use.
- Control solution results are accurate only between 59° – 104° F (15° – 40° C).
- Do *not* add water or other liquid to control solution.
- There are three levels of control solution available (Low, Normal, High).

For information on how to obtain control solution, contact Customer Care.

IMPORTANT: Control solution results should fall within the control solution range printed on the test strip vial. Do *not* use this range when testing blood glucose.

Out of Range Control Solution Results

- Repeat the test if control solution results are outside the range printed on the test strip vial.
- Stop using the meter if control solution results are consistently outside the range printed on the test strip vial. Contact Customer Care.

Causes of out-of-range results include:

- Expired or bad control solution
- Expired or bad test strip
- Testing error
- Watered-down control solution
- Meter malfunction
- Test strip deterioration
- Control solution test done outside 59° – 104° F (15° – 40° C)

Performing a Control Solution Test

1. Check control solution expiration date.


IMPORTANT: Do *not* use control solution past expiration date. Discard control solution 3 months after opening or on the expiration date printed on the bottle whichever comes first. (Example: open April 15, discard July 15; write the discard date on the side of the bottle.)

2. Remove test strip from vial.
3. Insert test strip to turn meter on.

Note: The meter turns off after 2 minutes of inactivity. Remove and reinsert the unused test strip to restart the meter.

If the meter does not turn on go to Troubleshooting section.

4. Confirm System Check Screen.

- This display appears when you turn on the meter. **Do not use the meter if the display check screen does not exactly match the example. Contact Customer Care.**
- See The FreeStyle Lite® Meter and How It Works section for more information.
- If the system check screen does not appear, go to Troubleshooting section.


5. Apply and **hold** control solution to 1 sample area of the test strip until:

- You see short lines moving clockwise on the display or
- You hear a beep

This indicates the test strip obtained enough control solution and the meter is processing. If the test does not start after applying the control solution, go to Troubleshooting section.

IMPORTANT: If a short line does not appear after 5 seconds, the sample may be too small. Do *not* add control solution to both sides of the strip. The test strips are for single use only. Do *not* reuse test strips. Discard used test strips.

(Example only)


6. View result.

The result appears on the display when the test is complete.

- ## 7. Compare the control solution result to the range printed on the test strip vial.
- The result should fall within the range.

Note: If you see an “Er1, Er2, Er3, or Er4” on the display, it is an error code. Go to Error Codes section.

(Example only)

- ## 8. Press and **hold** **C** for 2 seconds until appears on the display to mark the result as a control solution result. To unmark a control solution test result, press and **hold** **C** for 2 seconds. A value with the indicates you have marked this reading as a control solution test.

- ## 9. Remove the test strip to turn meter off.

Using Your Meter Memory

Purpose


The purpose of this section is to review your stored test results and 7-, 14-, and 30-day averages.

Viewing Blood Glucose Averages in Memory

Start with the meter off (no test strip inserted). Press **(m)**. Your 7-day average appears (control solution, LO, and HI readings *not* included).

Note: To exit meter memory at any time, press and hold **(m)**.

Example:

Memory Display

your average

n= number of results included in your average

- Press **(C)** to scroll through 14- and 30-day averages and all stored test results. The most recent result appears first.
- Press **(m)** to scroll backward.

Viewing Test Results in Memory


Meter memory stores the last 400 test results with time and date. Press  to scroll through the results. The most recent result appears first.

mem	on the display indicates that you are viewing stored test results.
davg	shows the average value of glucose readings from the last 7-, 14-, or 30-days.
	on the display indicates that the test result occurred when it was too hot or too cold to use the meter and may not be accurate. This result is included in 7-, 14-, and 30-day averages.
LO	a test result of under 20 mg/dL is included in memory as LO.
HI	a test result of over 500 mg/dL is included in memory as HI.
mem 	indicates a control solution test result.

Setting Reminder Alarms (optional)

Purpose

The purpose of this section is to set up to four reminder alarms to prompt you to test throughout the day.

Reminder Alarm

1. With the meter off, press and **hold** (m) for **3** seconds until SET appears and **0** is flashing.
2. Press (C). Reminder alarm 1 flashes.

Hour

Note: You are setting time for Reminder Alarm not clock time.

3. Set hour
 - Press (m) to set hour
 - Press (C) to display correct hour
 - Press (m) to save

Minutes

4. Set minutes
 - Press (C) to display correct minutes
 - Press (m) to save

5. Set alarm sound volume

If you prefer:	Press C until you see:	Then
No sound		• Press m to save
A low beep		
A louder beep		

6. Press **C** and repeat 2, 3, 4 and 5 to set remaining reminder alarms or

7. Press and **hold m** for 3 seconds until meter turns off.

Transferring Test Results To a Computer



PC


appears on the display when you insert a data cable in the data port. The computer must have a compatible data management system to complete the transfer. Remove this cable before performing a test.

Contact Customer Care to learn more about the data management system.


Replacing the Battery

The meter comes with a #2032, 3-volt, lithium battery installed. It provides power for about 500 tests.

- When  appears on the display, the battery is low.
- When  appears, the display light and test light dim, or fail to turn on.

IMPORTANT: When  appears, or when the display light and test light dim, replace the battery immediately.

- When you are ready to install a new battery, remove the old battery and install the new battery within 60 seconds to avoid losing time and date settings.
- If time and date settings are lost, both time 12:00 and date 00-00 will flash until reset. Reset time and date to get the correct 7-, 14-, and 30-day averages. If the time and date settings are lost, you may still perform an accurate blood glucose test.
- Memory log and previously stored user settings, such as date and time format, are not affected by removing the battery.

1. Press in the battery door clip and move the battery door away from the meter.
2. Push the battery toward the bottom of the meter to remove it.

3. Insert a new battery: (+) facing up.

4. To close, angle the battery door back onto the meter and snap into place.

Note: *The battery in this product should be removed and disposed in accordance with local regulations for separate collection of spent batteries.*

Cleaning the Meter

Avoid getting dirt, dust, blood, control solution, or liquid in the meter test strip and data ports.

Clean the outside of the meter using a damp cloth and:

- Mild detergent/soap and water, or
- 70% isopropyl alcohol, or
- A mixture of 1 part household bleach, 9 parts water

IMPORTANT:

Do *not* immerse the meter in water or other liquid.

Specifications

Assay method	Coulometric electrochemical sensor
Automatic shutoff	2 minutes of inactivity
Battery life	500 tests
Calibration	Plasma equivalent
Hematocrit	15% to 65%
Measurement units	mg/dL
Meter storage temperature	- 4° to 140° F (-20° to 60° C)
Memory	400 blood glucose and control solution tests with date and time
Operating relative humidity	5% to 90% (non-condensing)

Specifications (continued)

English

Operating temperature	40° to 104° F (4° to 40° C)
Power source	One CR 2032, 3V lithium battery, replaceable
Result range	20 to 500 mg/dL
Sample	Fresh capillary and venous whole blood samples
Sample size	0.3 microliter (300 nanoliters)
Size	1.57 in (w) x 2.9 in (l) x 0.65 in (d) 4.0 cm (w) x 7.4 cm (l) x 1.7 cm (d)
Weight	1.0 to 1.1 oz. (28.3 to 31.2 g) including battery

Definition of Symbols

	Consult instructions for use		For <i>in vitro</i> diagnostic medical device
	Temperature limitation		Catalog number
	Manufacturer		Recycle
	Batch code		Serial number
	Use by		Do not reuse
	The battery in this product should be removed and disposed in accordance with local regulations for separate collection of spent batteries.		

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©2010 Abbott
FreeStyle® and other brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions.
Abbott Diabetes Care Inc.