

# SELECT A BODYPOINT HIP BELT IN 6 STEPS

## 1. SELECT BELT TYPE

- Padded (pg. 5)
- Non-Padded (pg. 9)
- Hook-and-Loop-Compatible (pg. 8)

## 2. TWO OR FOUR-POINT



### TWO-POINT HIP BELT

Two points of attachment to the wheelchair. Suitable when the belt is supporting users with low tone or weakness.



### FOUR-POINT HIP BELT

Four points of attachment to the wheelchair. Primary straps position the belt like two-point hip belt and secondary straps anchor the belt in position to prevent it from riding up into the abdomen or twisting.


## 3. SELECT BUCKLE

OPTIONAL COVERS:



### Push-Button Buckle

Low profile, smooth, snag-free corners. Removable cover for changing release hole size and easy cleaning. Standard Access Cover (a) included. Reduced Access Cover (b) and Security Cover (c) sold separately.



### Plastic Side-Release

Lightweight, secure snap action.



### Rehab Latch™ Buckle

Lightweight, strong, low profile, easy to operate with limited hand function.



### Metal Push-Button

Strong, secure, child-sized buckle. Buckle Security Cover also available.

## 4. SELECT A PULL TYPE



### CENTER-PULL HIP BELT

**Operator:** User or caregiver

**Adjustment:** Tightens at the center buckle.

**Buckles:** Available in the full buckle range.

**Applications:** Clients with low muscle tone or weakness.



### REAR-PULL HIP BELT

**Operator:** Caregiver

**Adjustment:** Tightens from the rear of the pad pulling toward the center.

**Buckles:** Available in full buckle range

**Applications:** Clients with high tone. Counters pelvic rotation by allowing the caregiver to position one side of the pelvis, lock it in place and then position on the other side.



### DUAL-PULL HIP BELT

**Operator:** User or caregiver

**Adjustment:** Tightens from both sides of the center buckle.

**Buckles:** Side-release buckle only

**Applications:** Clients requiring greater adjustment to accommodate clothing and weight changes.

## 5. SELECT A SIZE

Sizes available from extra small to large.

See measurement section to the left for details.

## 6. SELECT A MOUNTING OPTION



**Cinch-Mount™**  
Time-saving, precise, easy adjustment


**Flat-Mount™**  
Versatile, lowest profile



### FrameSaver™ Clamp (optional)


Create a mounting point without drilling. Fits Ø7/8"–Ø1" (22-25mm) tube.

## MEASURING FOR A HIP BELT



Measure hip width across the greater trochanters with the person seated. If the hip width falls between sizes on our chart, consider other factors such as growth, weight changes, and clothing.

<b>XS</b>	5–9" (13–23cm)
<b>S</b>	7–11" (18–28cm)
<b>M</b>	9–15" (23–38cm)
<b>L</b>	13–19" (33–48cm)



	A	B	C	D
<b>XS</b>	1" (25mm)	1-3/4" (44mm)	6" (15cm)	50" (127cm)
<b>S</b>	1-1/2" (38mm)	2-1/4" (57mm)	7" (18cm)	55" (140cm)
<b>M</b>	1-1/2" (38mm)	2-1/2" (64mm)	9" (23cm)	60" (152cm)
<b>L</b>	2" (50mm)	3" (76mm)	11" (28cm)	65" (165cm)