OneWaySlide

SystemRoMedic™

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2074
2075
2076
2077

SWL: 150 kg/330 lbs
200kg/440 lbs

OneWaySlide is used for repositioning farther back and sideways in a chair/wheelchair. OneWaySlide is designed to glide easily backward, but not forward. To make manoeuvring easier for caregivers, OneWaySlide features handles on both sides and at the back. It is placed so that the user remains sitting on it. OneWaySlide is available in several sizes.


Check Safety

Visual inspection


Check the condition and function of the assistive device regularly. Always inspect the product after laundering. Check to ensure that seams and material are free from damage. Check to ensure that fabric is not worn or faded. If there are signs of wear, the product must be discarded.

Mechanical load

Apply pressure to buckles and handles and verify that they will not open or break under heavy loads.

Friction test

Low friction - Ensure that the surface is smooth and slides effortlessly. Compare with a new product. →

High friction - Ensure that the surface is rough and prevents the transfer aid from moving. Compare with a new product. ←


Always read the user manual

Always read the user manuals for all assistive devices used during a transfer.


Keep the user manual where it is accessible to users of the product.

Always make sure that you have the right version of the user manual.

Placement of OneWaySlide



Place OneWaySlide so that the open section is facing the back of the wheelchair. Place the upper section so that it sticks up a bit on the backrest of the wheelchair. This makes it easier for the caregiver to reach the handles.


Tie the OneWaySlide to the rear frame tubing of the wheelchair.

Alternative grips


The caregiver can either simply grasp the handle, or first loop it around the wrist. Remember to work with knees bent, back and arms straight, and shift the weight.


Using OneWaySlide


Positioning is easier if the patient can lean forward. If the user is unable to lean forward, one caregiver can assist from the front while supporting the patient's balance.

If the user can lean forward and has good balance while seated, only one caregiver is required. Grasp the appropriate handles and pull the user back by drawing the OneWaySlide upward along the backrest of the wheelchair.

If the user is heavy or if the chair has a high back, two caregivers can assist, one working on either side. Coordinate your efforts.

A caregiver can easily help make small positioning changes sideways using the handles on each side of the OneWaySlide

! Raising the user's feet from the floor will relieve some of the weight of the thighs from the seat, thereby making repositioning easier.

Material

Polyester, polyamide, cotton

Care of the product

Read the product label.

Do not use rinse/ fabric softener and avoid tumbler-drying for optimal service life.

CE Medical Device Class I

Medical Device Class I. The product complies with the requirements of the Medical Devices Directive 93/42/EEC.

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Simple solutions for great results

SystemRoMedic™ is the name of Handicare's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic™ is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic™ is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic™ offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

Always make sure that you have the correct version of the manual

For questions about the product and its use

Please contact your local Handicare and SystemRoMedic™ representative.

Handicare offers solutions and support to increase the independence of disabled or elderly people as well as to improve the convenience of those who are caring for them.

The Handicare Group is one of the leading healthcare companies in Europe with its own manufacturing organizations and sales companies in Norway, Sweden, Denmark, Germany, the Netherlands, Great Britain, France, China, Canada and the USA. Handicare's products are also distributed by partners in more than 40 countries worldwide. Our wide range of high-quality products includes various manual and power wheelchairs, seating systems, scooters, a complete easy transfer system and other patient handling aids, stairlifts, car adaptations, rise and recline chairs, and bathing and toileting products.

Handicare AB
Maskinvägen 17
SE-972 54 Luleå, SWEDEN