









**Patented Breath Sensing** 

CTOX-MN02

The SmartDose® Mini Auto-Adjusting Oxygen Conserver senses patient breathing patterns, and adjusts the preset rest setting higher by up to 2 settings as needed to best match dose to need.

Auto-adjusting conservers are the future of oxygen therapy, optimizing both patient saturations and oxygen conservation. Satisfy your patients' and referrals' wishes while decreasing your delivery costs, all with one unit.

# **SmartDose® Mini Auto-Adjusting Oxygen Conserver**

Item #	Product Description		UOM
CTOX-MN02	SmartDose® Mini Auto-Adjusting Oxygen Conserver		1/ea
		•	Allows delivered dose to increase up to 2 settings during patient exertion up to a maximum pulse dose
The state of the s			setting of 7. Can use 2 AA Alkaline or 2 AA rechargeable batteries
A chiovos ma	vimum conconvation officions v by		for a duration of about 1 year for most nationts *

- Achieves maximum conservation efficiency by automatically adjusting the dose down when less oxygen is needed.
- Delivers a quiet, gentle dose of oxygen using an efficient and compact design.
- for a duration of about 1 year for most patients.3
- Operating pressure range: 500 3000 PSI.
- Weighs just 15.5 ounces with batteries.
- Auto shut-off feature conserves battery life.
- Fits easily in virtually any carry bag.

\*Average 4 hours per day usage.

accessories		UOM
СТВТ-04СН	4 Rechargeable Batteries and Battery Charger for SmartDose® Mini and Liquid Oxygen Units	1/ea
SEAL100	Replacement Seal-Washer for Oxygen Regulators; Made from Brass and Viton	10/bag



# **DURATION CHARTS**

#### **CHAD® Bonsai® Velocity**

PNEUMATIC OXYGEN CONSERVER

				Sett	ting				CF
Cylinder Type	Cylinder Volume	1	2	3	4	5	6	CF	
			Estimated Cylinder Duration in Hours (Based on 20 breaths per minute)						4 lpm
M4 (A)	113 Liters	7.2	4.3	2.9	2.4	2.0	1.8	0.9	0.5
M6 (B)	164 Liters	10.5	6.2	4.3	3.5	3.0	2.6	1.4	0.7
ML6	171 Liters	11.0	6.5	4.5	3.7	3.1	2.7	1.4	0.7
M9 (C)	246 Liters	15.8	9.3	6.4	5.3	4.5	3.9	2.1	1.0
D	425 Liters	27.2	16.1	11.1	9.1	7.7	6.8	3.5	1.8
E	680 Liters	43.6	25.8	17.7	14.5	12.3	10.9	5.7	2.8

# **CHAD®** Evolution™

# **CHAD®** Evolution<sup>™</sup> Motion

ELECTRONIC OXYGEN CONSERVER

AUTO-ADJUSTING OXYGEN CONSERVER

ELECTROTTIC C	ACTO ADSOSTINO OXTOLINO								POTTOLICALI	
		Setting								
Cylinder Type	Cylinder Volume	1	2	3	4	5	6	7	CF	
		1	Estimated Cylinder Duration in Hours (Based on 20 breaths per minute.  Durations shown are for base setting. Auto-adjusting feature may shorten times)							
M4 (A)	113 Liters	9.4	4.7	3.1	2.4	1.9	1.6	0.7	0.9	
M6 (B)	164 Liters	13.7	6.8	4.6	3.4	2.7	2.3	1.9	1.4	
ML6	171 Liters	14.3	7.1	4.8	3.6	2.9	2.4	2.0	1.4	
M9 (C)	246 Liters	20.5	10.3	6.8	5.1	4.1	3.4	2.9	2.1	
D	425 Liters	35.4	17.7	11.8	8.9	7.1	5.9	5.0	3.5	
E	680 Liters	56.7	28.3	18.9	14.2	11.3	9.4	8.1	5.7	

#### **SmartDose® Mini**

AUTO-ADJUSTING OXYGEN CONSERVER

		Setting							
Cylinder Type	Cylinder Volume	1	2	3	4	5	CF		
			Estimated Cylinder Duration in Hours (Based on 20 breaths per minute.  Durations shown are for base setting. Auto-adjusting feature may shorten times)						
M4 (A)	113 Liters	5.9	2.9	2.0	1.5	1.2	0.9		
M6 (B)	164 Liters	8.5	4.3	2.8	2.1	1.7	1.4		
ML6	171 Liters	8.9	4.5	3.0	2.2	1.8	1.4		
M9 (C)	146 Liters	7.6	3.8	2.5	1.9	1.5	2.1		
D	425 Liters	22.1	11.1	7.4	5.5	4.4	3.5		
E	680 Liters	35.4	17.7	11.8	8.9	7.1	5.7		

#### **SmartDose® LOX**

LIQUID OXYGEN PORTABLE CONSERVER SYSTEM

		Setting									
Cylinder Type	Cylinder Volume	1	2	3	4	5	CF				
			Estimated LOX Duration in Hours (Based on 20 breaths per minute.  Durations shown are for base setting. Auto-adjusting feature may shorten times)								
0.3 LOX	285 Liters	14.9	14.9 7.4 5.0 3.7 3.0								
0.6 LOX	540 Liters	28.1	14.1	9.4	7.0	5.6	4.5				