Care and Application

for Compression Stockings and Pantyhose

Freedom in Motion™

Thank You!

Thank you for choosing Juzo® compression stockings. Manufactured with care and detail, Juzo® has been an international leader in providing high quality medical compression garments since 1912. Most Juzo® stockings feature high quality Lycra® covered with specially selected textile fibers. These advanced knitting techniques and superior raw materials enable Juzo® to provide garments that offer optimal wearing comfort, air permeability, and largely resistance to sun, ointments, lotions and perspiration.

Care and Maintenance

To maintain your Juzo® stockings, we recommend the following:

- Hand or machine wash warm daily to maximize therapeutic effectiveness and prolong the life of your Juzo® garment.
- When machine washing and drying, use the lowest permanent press settings. A garment bag is strongly recommended and available in the *Juzo® Compression Therapy Care accessory package*.
- Wash dark colors separately from other clothing.
- Use Juzo® washing detergent or another mild soap.
- DO NOT USE: bleach, chlorine, fabric softeners, or other laundry additives. These may damage your garment.
- Please refer to the washing symbols on the label sewn into your Juzo stocking.
- Store your stockings in a dry place at room temperature.

Return Policy

Should it become necessary to return your garment for evaluation, for hygienic reasons Juzo® requires that the garment be laundered prior to returning to the original place of purchase.

Fabric

Details are provided on label sewn inside your compression stocking.

Instructions for Use

It is recommended that you consult with your medical professional and fitter for proper guidelines on the use and application of this product. As therapy progresses or new symptoms develop, it may be necessary to take new measurements for a new compression stocking.

It is also recommended that you put your compression stockings on first thing in the morning every day. Rough fingernails, jagged rings etc. can damage fabric; use a Juzo latex free donning glove to protect garment.

Indications

Primary & secondary varices, pregnancy varices and or edema, venous & chronic venous insufficiency, post superficial thrombophlebitis, prophylaxis of thrombosis & the post-thrombotic conditions, post-thrombotic syndrome, postoperative & post-traumatic edema, lymphatic edema & lipedema, elephantiasis, angiodysplasia, pre & post vein procedure, arthrogenous congestive syndrome, paresis & partial paresis of extremities, prophylaxis of venous disorders, heavy & fatigued legs, atrophy blanche, dermatosclerosis, post ulcer, severe tendency to edema & treatment of burn scars.

Contraindications*

Absolute: peripheral arterial disease, decompensated heart diseases, septic phlebitis, phlegmasia coerulea dolens.








Relative: concomitant dermatoses, intolerance to compression stocking fabric, sensory disturbance of the limb, advanced peripheral neuropathy, primary chronic arthritis.

Cautions

No side effects are currently known. However if a negative result should occur (i.e. skin irritations) from the use of our product, please consult with a medical professional. The manufacturer is not liable for damage / injuries which have been caused by improper handling or misuse.

Julius Zorn Inc assumes no liability as a result of any contraindicated use of this product.

*Contraindications according to the German Phlebology Society.

Useful tips to help with your compression stocking:

For knee-highs, thigh-highs and pantyhose:

- 1 Turn the upper half of the stocking inside out to the heel by gently sliding the stocking down over the foot portion.
- 2 Slip your foot into the foot portion.
- 3 Gently ease the remaining part of the stocking over the heel.
- 4 Use the palms of your hands to massage the stocking upwards towards the knee. *NOTE: Never pull at the upper edge of the stocking.*
- 5 Follow the same procedure for the opposite leg, taking care that the stitch rows always run vertically.
- 6 **For thigh-highs and pantyhose**, position both stockings, as described above, until just below the knee
- 7 Gently ease the stocking onto the thigh. Work the stocking upwards into place using the palms of your hands until the stocking reaches the groin area.
NOTE: If putting on a thigh-high stocking with silicone border, never pull the stocking up from the silicone border.
- 8 **For pantyhose**, ease the body part of the panty hose to your waist. The seams should run vertically up the front of the garment.

Always use Juzo latex free donning gloves when applying compression garments.

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Accessories

Juzo® Detergent

For laundering compression garments and other fine knits.

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Juzo® Slippie Gator

Juzo® Slippie Gator was designed to be a lightweight, easy to use, donning and doffing device for both open and closed toe compression stockings/socks. When used with the Juzo® Easy Pad, the Slippie Gator provides nearly hands-free donning and doffing, minimizing bending at the waist.

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Juzo® Adhesive Lotion

Special roll-on lotion to keep your garment in place. Use as an alternative to other forms of attachments.

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Juzo® Latex Free Gloves

Provide excellent protection against rings and fingernails, and help you to correctly position your garment.

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Alps Fitting Lotion with Skin Conditioner

100% Silicone Fitting Lotion which helps prevent snags and facilitates the donning of compression garments.

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