# **AUTOMATIC WRIST BLOOD PRESSURE UNIT**

**(ELECTRONIC SPHGMOMANOMETER)**  
**OPERATION GUIDE**

Read all of the information in the operation guide and any other literature in the box before operating the unit.

**REF** MDS3003    **REF** MDS3003P

# CONTENTS

<b>Important Information</b>	3	7. Displaying Stored Results	17
<b>Contents And Display Indicators</b>	4	8. Deleting Measurements From The Memory	19
<b>Intended Use</b>	5	9. Assessing High Blood Pressure For Adults	19
<b>Contraindication</b>	5	Troubleshooting (1)	21
<b>Product Description</b>	6	Troubleshooting (2)	22
<b>Specifications</b>	7	<b>Maintenance</b>	24
<b>Notice</b>	8	<b>Explanation Of Symbols On Unit</b>	25
<b>Setup And Operating Procedures</b>	12	<b>Warranty Information</b>	26
1. Battery Loading	12	<b>Service Center</b>	26
2. Clock And Date Adjustment	12	<b>Electromagnetic Compatibility Information</b>	28
3. Connecting The Cuff To The Monitor	14		
4. Applying The Cuff	14		
5. Body Posture During Measurement	15		
6. Taking Your Blood Pressure Reading	16		

## IMPORTANT INFORMATION

### **NORMAL BLOOD PRESSURE FLUCTUATION**


All physical activity, excitement, stress, eating, drinking, smoking, body posture and many other activities or factors (including taking a blood pressure measurement) will influence blood pressure value. Because of this, it is mostly unusual to obtain identical multiple blood pressure readings.

Blood pressure fluctuates continually day and night. The highest value usually appears in the daytime and the lowest one usually at midnight. Typically, the value begins to increase at around 3:00 AM, and reaches the highest level in the daytime while most people are awake and active.

Considering the above information, it is recommended that you measure your blood pressure at approximately the same time each day.

Please always relax a minimum of 1 to 15 minutes between measurements to allow the blood circulation in your arm to recover. It is rare that you obtain identical blood pressure readings each time.

## CONTENTS AND DISPLAY INDICATORS

- A. Memory Indicator
- B. Low Battery Indicator
- C. Blood Pressure Level Classification Indicator
- D. Irregular Heartbeat Symbol
- E. Ready to Inflate Indicator
- F. Date/Time Display (Alternating)
- G. Systolic Pressure
- H. Diastolic Pressure/Pulse Rate Display (Alternating)

- A. WHO Color Bar
- B. LCD Display
- C. Cuff
- D. MEM Button
- E. START Button

## INTENDED USE

Fully Automatic Electronic Blood Pressure Monitor is for use by medical professionals or at home and is a non-invasive blood pressure measurement system intended to measure the diastolic and systolic blood pressures and pulse rate of an adult individual by using a non-invasive technique in which an inflatable cuff is wrapped around the wrist.

## CONTRAINDICATION

⚠ It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.

## PRODUCT DESCRIPTION

Based on Oscillometric methodology and silicon integrated pressure sensor, blood pressure and pulse rate can be measured automatically and non-invasively. The LCD display will show blood pressure and pulse rate. The most recent 1 X 60 measurements can be stored in the memory with date and time stamp. The monitor can also show the average reading of the last three measurements. The Electronic Sphygmomanometer corresponds to the below standards: IEC 60601-1:2005(Medical electrical equipment-Part1:General requirements for safety) , IEC60601-1-2:2007(Medical electrical equipment-Part1:General requirements for safety; Collateral Standard-Electromagnetic compatibility-Requirements and tests), EN 1060-1:1995+A1:2002+A2:2009 (Non-invasive sphygmomanometers - Part 1: General requirements), EN 1060-3:1997+A1:2005+A2:2009 (Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems). IEC 80601-2-30: 2009+Cor.2010(Medical electrical equipment -Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers)

## SPECIFICATIONS

1. Product name: Blood Pressure Monitor
2. Model: MDS3003/MDS3003P
3. Classification: Internally powered, Type BF applied part, IPX0, No AP or APG, Continuous operation
4. Machine size: 83mm x 64mm x 28mm (3 1/4" x 2 17/32" x 1 3/32")
5. Cuff circumference: 14cm - 19.5cm (5 1/2" - 7 11/16"),
6. Weight: approx. 110g (3 7/8oz.) (excludes batteries)
7. Measuring method: oscillometric method, automatic air inflation and measurement
8. Memory volume: 60 times with time and date stamp
9. Power source: 2 × 1.5V  $\Rightarrow$  SIZE AAA batteries
10. Cuff pressure: 0-300mmHg
11. Range of measurement:
  - Systolic: 60-260mmHg
  - Diastolic: 40-199mmHg
12. Measuring accuracy:  $\pm 3$ mmHg
13. Pulse rate range:  $\pm 5\%$

14. Environmental temperature for operation: 5°C~40°C (41°F~104°F)
15. Environmental humidity for operation: ≤90% RH
16. Environmental temperature for storage and transport: -20°C~55°C (-4°F~131°F)
17. Environmental humidity for storage and transport: ≤90% RH
18. Environmental pressure: 80KPa-105KPa
19. Battery life: Approx. 270 measurements.
20. List of all components belonging to the pressure measuring system, including accessories: Pump, Valve, LCD, Cuff, Sensor

## NOTICE

1. Keep still, calm and rest for 5 minutes before blood pressure measurement.
2. The cuff should be placed at the same level as your heart.
3. During measurement, neither speak nor move your body and arm.
4. Measure on the same wrist each measurement.
5. Consult with your physician for any of the instances below:  
—This device does not measure circulation or blood flow.

- The application of the cuff is over a wound or inflammation.
  - The application limb shares intravascular access or therapy, or an arteriovenous (A-V) shunt.
  - The application arm is on the same side as a mastectomy.
  - The application limb shares monitoring equipment.
6. Please always relax a minimum of 1 to 15 minutes between measurements to allow the blood circulation in your arm to recover. Prolonged over-inflation (cuff pressure exceeding 300 mmHg or maintained above 15 mmHg for longer than 3 minutes) of the bladder may cause ecchymoma of your arm.
  7. Remove the batteries if the monitor will not be used for a month or more to avoid damage caused by battery leakage.
  8. ⚠️ This Electronic Sphygmomanometer is designed for adults and should never be used on infants or young children. Consult your physician or other health care professionals before use on older children.
  9. Blood pressure measurements determined by this monitor are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard Institute, electronic or automated sphygmomanometers.

10. Information regarding potential electromagnetic or other interference between the blood pressure monitor and other devices together with advice regarding avoidance of such interference please see part ELECTROMAGNETIC COMPATIBILITY INFORMATION.
11. If Irregular Heartbeat (IHB) brought by common arrhythmias is detected in the procedure of blood pressure measurement, a signal of '(♥)' will be displayed. Under this condition, the Electronic Sphygmomanometer can keep function, but the results may not be accurate; it's suggested that you consult with your physician for accurate assessment.

There are 2 conditions under which the signal of IHB will be displayed:

- 1) The coefficient of variation (CV) of pulse period  $>25\%$ .
  - 2) The difference of adjacent pulse period  $\geq 0.14s$ , and the number of such pulse takes more than 53 percent of the total number of pulse.
12. Please do not use the cuff other than supplied by the manufacturer; otherwise it may bring biocompatible hazard and might result in a measurement error.
  13. ⚠ The monitor might not meet its performance specifications or cause safety hazard if stored or used outside the specified temperature and humidity ranges in specifications.

14. ⚠ Please do not share the cuff with other infective person to avoid cross-infection.
15. Please note that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
16. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
  - Reorient or relocate the receiving antenna.
  - Increase the separation between the equipment and receiver.
  - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
  - Consult the dealer or an experienced radio/TV technician for help.

## SETUP AND OPERATING PROCEDURES

### 1. BATTERY LOADING


- a. Open battery cover at the back of the monitor.
  - b. Load 2 “AAA” size batteries. Please pay attention to polarity.
  - c. Close the battery cover.
- ⚠ When LCD shows battery symbol , replace all batteries with new ones.
- ⚠ Rechargeable batteries are not suitable for this monitor.
- ⚠ Remove the batteries if the monitor will not be used for a month or more to avoid damage or battery leakage.
- ~~⚠~~ The monitor, the batteries and the cuff, must be disposed of according to local regulations at the end of their usage.

### 2. CLOCK AND DATE ADJUSTMENT


- a. Once you install the battery or turn off the monitor, it will enter Clock Mode, and LCD will display time and date by turns. See picture 2 and 2.1.
- b. While the monitor is in Clock Mode, press both the “START” and “MEM” button simultaneously, a beep is heard and the month will blink at first. See picture 2-2. Press the button “START” repeatedly, the day, hour and minute will blink turn.

While the number is blinking, press the button “MEM” to increase the number. Keep on pressing the button “MEM”, the number will increase fast.


- c. You can turn off the monitor by pressing “START” button when the minute is blinking, then the time and date is confirmed.
- d. The monitor will turn off automatically after 1 minute of no operation, with the time and date unchanged.
- e. Once you change the batteries, you should readjust the time and date.

Picture 2


Picture 2.1

Picture 2.2


### 3. CONNECTING THE CUFF TO THE MONITOR

The cuff is attached to the monitor when it is packaged, Should the cuff become unattached, align the two plugs and four brackets of the cuff with the plug sockets and bracket sockets of the monitor and press the cuff to the monitor until the plugs and brackets are securely attached.

### 4. APPLYING THE CUFF

- Place the cuff around a bare arm 1-2 cm above the wrist joint on the palm side of the wrist.
- While seated, place the arm with the cuffed wrist in front of your body on a desk or table with the palm up. If the cuff is correctly placed, you can read the LCD display.
- The cuff must be neither too tight nor too loose.

## Note:

- Please refer to the cuff circumference range in “SPECIFICATIONS” to make sure that the appropriate cuff is used.
- Measuring on the same wrist each time.
- Do not move your arm, body, or the monitor during measurement.
- Keep still, calm for 5 minutes before blood pressure measurement.
- Please keep the cuff clean. If the cuff becomes dirty, remove it from the monitor and clean it by hand in a mild detergent, then rinse it thoroughly in cold water. Never dry the cuff in clothes dryer or iron it. Clean the cuff after the usage of every 200 times is recommended.
- Do not place the cuff around your wrist if the wrist has any inflammation, acute diseases, infections skin wounds.

## 5. BODY POSTURE DURING MEASUREMENT

### Sitting Measurement

- a. Be seated with your feet flat on the floor.
- b. Place palm upside in front of you on a flat surface such as desk or table, with your elbow resting on a chair or table.
- c. The cuff should be at the same level as your heart.

- d. Sit comfortably and don't cross your legs.
- e. Placing your arm on the monitor's plastic case as shown in the picture to aid in maintaining a correct and still arm position.

## 6. TAKING YOUR BLOOD PRESSURE READING

- a. After applying the cuff with your body is in a comfortable position, press the "START" button. A beep is heard and all display characters are shown for self-test. See picture 6. Please contact the service center if segment is missing.

Picture 6

Picture 6.1

Picture 6.2

Picture 6.3

- b. Then the monitor starts to seek zero pressure. See picture 6-1.
- c. The monitor inflates the cuff until sufficient pressure has built up for a measurement. Then the monitor slowly releases air from the cuff and carries out the measurement. Finally the blood pressure and pulse rate will be calculated and displayed on the LCD screen separately. Irregular heartbeat symbol (if any) will blink. See picture 6.2 and

- 6.3. The result will be automatically stored in the memory bank.
- d. After measurement, the monitor will turn off automatically after 1 minute of no operation. Alternatively, you can press the “START” button to turn off the monitor manually.
  - e. During measurement, you can press the “START” button to turn off the monitor manually.


**Note: Please consult a healthcare professional for interpretation of pressure measurements.**

## **7. DISPLAYING STORED RESULTS**

- a. After measurement, you can review the results in the memory bank by pressing the “MEM” button. Alternatively, you can press “MEM” button in Clock Mode to display the stored results. If it no result stored, LCD will show dashes as picture 7, while press the button “MEM” or “START”, machine will turn off. If there are results in the memory bank, the LCD will display the amount of the results in the memory bank. See picture 7.1.
- b. Then, the most recent result will be displayed with date and time stamp. See picture 7.2. The blood pressure and pulse rate will be shown separately. Irregular heartbeat symbol (if any) will blink. See picture 7.3 and 7.4. Press “MEM”

pressing the MEM button displays the respective results measured previously.

- c. When displaying the stored results, the monitor will turn off automatically after 1 minute of no operation. You can also press the button "START" to turn off the monitor manually.

Picture 7

Picture 7.1

Picture 7.2

Picture 7.3


Picture 7.4

Picture 7.5

## 8. DELETING MEASUREMENTS FROM THE MEMORY

When any result is displayed, keep on pressing button “MEM” for three seconds, all results in the memory bank will be deleted after three “beeps”. LCD will show dashes. See picture 8, press the button “MEM” or “START”, the monitor will turn off.


Picture 8

## 9. ASSESSING HIGH BLOOD PRESSURE FOR ADULTS

The following guideline for assessing high blood pressure (without regard to age or gender) has been established by the World Health Organization (WHO). Please note that other factors (e.g. diabetes, obesity, smoking, etc.) need to be taken into consideration. Consult with your physician for accurate assessment, and never change your treatment by yourself.

## Classification of blood pressure for adults

BLOOD PRESSURE CLASSIFICATION	SBP mmHg	DBP mmHg	COLOR INDICATOR
Optimal	<120	<80	GREEN
Normal	120-129	80-84	GREEN
High-Normal	130-139	85-89	GREEN
Grade 1 Hypertension	140-159	90-99	YELLOW
Grade 2 Hypertension	160-179	100-109	ORANGE
Grade 3 Hypertension	$\geq 180$	$\geq 110$	RED

*WHO/ISH Definitions and Classification of Blood Pressure Levels*

**Note: It is not intended to provide a basis of any type of rush toward emergency conditions/diagnosis based on the color scheme and that the color scheme is meant only to differentiate between the different levels of blood pressure.**

## TROUBLESHOOTING 1

<b>PROBLEM</b>	<b>POSSIBLE CAUSE</b>	<b>SOLUTION</b>
LCD Display shows abnormal result	The cuff position was not correct or it was not properly tightened	Apply the cuff correctly and try again
	Body posture was not correct during testing	Review the "BODY POSTURE DURING MEASUREMENT" sections of the instructions and re-test
	Speaking, arm or body movement, angry, excited or nervous during testing	Re-test when calm and without speaking or moving during the test
	Irregular heartbeat (arrhythmia)	It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer

## TROUBLESHOOTING 2

<b>PROBLEM</b>	<b>POSSIBLE CAUSE</b>	<b>SOLUTION</b>
LCD shows low battery symbol 	Low Battery	Change the batteries
LCD shows "Er 0"	Pressure system is unstable before measurement	Don't move and try again
LCD shows "Er 1"	Fail to detect diastolic pressure	
LCD shows "Er 2"	Fail to detect systolic pressure	
LCD shows "Er 3"	Pneumatic system blocked or cuff is too tight during inflation	Apply the cuff correctly and try again
LCD shows "Er 4"	Pneumatic system leakage or cuff is too loose during inflation	

## TROUBLESHOOTING 2.2

PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD shows "Er 5"	Cuff pressure above 300mmHg	Measure again after five minutes. If the monitor is still abnormal, please contact the local distributor or the factory
LCD shows "Er 6"	More than 3 minutes with cuff pressure above 15 mmHg	
LCD shows "Er 7"	EEPROM accessing error	
LCD shows "Er 8"	Device parameter checking error	
LCD shows "Er 9"	MCU self-verify error	
LCD shows "Er A"	Pressure sensor parameter error	
No response when you press button or load battery	Incorrect operation or strong electromagnetic interference	Take out batteries for five minutes, and then reinstall all batteries.

## MAINTENANCE

1. Do not drop this monitor or subject it to strong impact.
2. Avoid high temperature and sunlight. Do not immerse the monitor in water as this will result in damage to the monitor.
3. If this monitor is stored near freezing, allow it to acclimate to room temperature before use.
4. Do not attempt to disassemble this monitor.
5. If you do not use the monitor for a long time, please remove the batteries.
6. It is recommended the performance should be checked every 2 years or after repair. Please contact the service center.
7. If the monitor becomes dirty, please clean it with a soft dry cloth. Do not use any abrasive or volatile cleaners.
8. No component in the monitor can be maintained by the user. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment that are designated as reparable can be supplied.
9. The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years.
10. It is recommended that the cuff is disinfected 2 times per week. Wipe the inner side of the cuff that makes direct contact with the skin with a soft cloth and Ethyl alcohol (75-90%). Let the cuff air dry before next use.
11. The cuff can maintain safety and performance characteristics for an average of 1,000 measurements.

## EXPLANATION OF SYMBOLS ON UNIT

Symbol for “THE OPERATING GUIDE MUST BE READ”

Symbol for “WARNING”

Symbol for “Cuff Type: BF Applied Part”

Symbol for “ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local Authority or retailer for recycling advice”.

### **IP22**

The first characteristic numeral symbol for “Degrees of protection against access to hazardous parts and against solid foreign objects”. The second characteristic numeral symbol for “Degrees of protection against ingress of water”. The bigger the number is, the higher the level of protection is.

## WARRANTY INFORMATION

The Medline automatic electronic blood pressure monitor MDS3003/MDS3003P is under warranty for 2 years from the date of purchase.

This warranty includes the instrument/cuff. The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the instrument. The warranty is only valid upon presentation of proof of purchase.

## SERVICE CENTER

Manufactured for

Medline Industries, Inc., Three Lakes Drive, Northfield, IL 60093 USA.

Made in China. Tel: 1-800-MEDLINE RC20AXA

# ELECTROMAGNETIC COMPATIBILITY INFORMATION

## For all ME EQUIPMENT and ME SYSTEMS

### Guidance and manufacturer's declaration - electromagnetic emissions

The [MDS3003]/[MDS3003P] is intended for use in the electromagnetic environment specified below. The customer or the user of the [MDS3003]/[MDS3003P] should assure that it is used in such an environment.

<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The [MDS3003]/[MDS3003P] uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions: CISPR 11	Class B	The [MDS3003]/[MDS3003P] is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions: IEC 61000-3-2	No Applicable	
Voltage fluctuations/ flicker emissions: IEC 61000-3-3	No Applicable	

## Guidance and manufacturer's declaration - electromagnetic immunity

The [MDS3003]/[MDS3003P] is intended for use in the electromagnetic environment specified below. The customer or the user of the [MDS3003]/[MDS3003P] should assure that it is used in such an environment.

<b>IMMUNITY test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 6$ kV contact $\pm 8$ kV air	$\pm 6$ kV contact $\pm 8$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power Frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration - electromagnetic immunity

The [MDS3003]/[MDS3003P] is intended for use in the electromagnetic environment specified below. The customer or the user of the [MDS3003]/[MDS3003P] should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the [MDS3003]/[MDS3003P], including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance:</b></p> $d = 1.2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$

Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p><b>a</b> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the [MDS3003]/[MDS3003P] is used exceeds the applicable RF compliance level above, the [MDS3003]/[MDS3003P] should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the [MDS3003]/[MDS3003P].</p> <p><b>b</b> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

**Recommended separation distances between  
portable and mobile RF communications equipment and the [MDS3003]/[MDS3003P]**

The [MDS3003]/[MDS3003P] is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the [MDS3003]/[MDS3003P] can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the [MDS3003]/[MDS3003P] as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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