XPD 12

User manual

EN

User Manual

English

Contents

2	Foreword	14	Indications of use
3	General information	15	Unique benefits
4	Overview	16	Regular application principles
5	Charging the XPD 12 battery	17	Regular application methods
6	Features and mode description	17	Pain in the trunk
7	Directions for use	18	Pain in the joints and limbs
8	Electrode guidelines	19	Electrode placement variations
9	Electrode recommendations	20	Troubleshooting
9	Maintenance and care	22	Contraindications
10	Environmental specifications	23	Warnings
10	Technical specifications	25	Precautions
11	Pulse parameters	25	Adverse reactions
12	Mode description, TENS	25	Location of stimulation
13	Mode description, EMS	26	Conformity to safety standards
		27	Warranty information

-1-

Foreword

Before operation, please read carefully this user's manual and be clear about the instructions.

Do you suffer from **PAIN, MUSCLE TENSION** or **STRESS**?

The HiDew XPD 12 can quickly relax tense muscles and relieve pain. It is modern day home therapy device using the latest bio-electrical technology.

The HiDew XPD 12 is designed to provide endless amount of relaxing, deep, soothing sensations in the areas you need relief. Relax, feel good and rejuvenate your body and soul. Don't live with unhealthy physical and mental tension anymore.

Enjoy the XPD 12 and live life at fullest, whether you're at home, at work or on vacation. You'll love how good it feels!

Please read the instructions carefully before using this instrument.

The XPD 12 is a double-channel battery operated muscle stimulation system that helps to relieve minor muscular aches and pains. It comprises of two main components: an electronic stimulatory module which generates the required stimulation signals and 2 sets of pre-impregnated ECG electrode pads with lead wires.

To be used for temporary pain relief associated with sore and aching muscles in the lower back, upper extremities (arms), lower extremities (legs) and joints due to strain from exercise or normal household work activities. It is intended to be used to stimulate healthy muscles in order to improve and facilitate muscle performance.

-2-

General information

Included in this package:

- XPD 12 unit
- 1 set of LARGE single sided adhesive electrode pads (2.4 inches - round shape)
- 1 set of EXTRA-LARGE single sided adhesive electrode pads (3 1/2 inch - rectangular shape)
- 2 Electrode Wires
- USB Cable
- AC Adapter
- User Manual
- Pad Holder

Structure:

- | | |
|--------------------|--------------------------|
| 1. Control unit | 7. Output |
| 2. USB cable | 8. LCD screen |
| 3. AC adapter | 9. Power switch |
| 4. Electrode wires | 10. USB port |
| 5. Electrode pads | 11. A/B Channel selector |
| 6. Pad holder | |

*In case of any discrepancy between the physical and picture, refer to the actual product.

-3-

Overview: control unit and accessories

USB cable: use to charge the XPD 12 battery

AC adapter: use with USB cable to charge battery

Electrodes: use to attach to your bare skin to perform stimulation

Electrode wires: use to make the connection between the device and the electrodes

Pad holder: use to store electrode pads and wires after use

-4-

Charging the XPD 12 battery

*Via a standard wall outlet

1. Connect the AC adapter to any standard wall outlet.

2. Connect the small end of the USB cable to the unit and the bigger end to the AC adapter.

3. The unit is finished charging when the battery icon indicates full.

*Via a computer USB port

1. Connect the small end of the USB cable to the unit and bigger end to the USB port on a computer (computer must be turned on).

2. The unit is finished charging when the battery icon indicates full.

-5-

*Notes:

1. The lithium battery can be recharged through both AC adapter and USB cable.
2. While charging, the screen in ON position indicates the charge, while in OFF position the screen is blank. However the battery will charge in either one of the power positions (ON or OFF).
3. The battery should be charged for 8-10 hours before first use.

Features

FDA approved Class II Medical Device

- A large, easy to read LCD display
- Backlight for increased readability
- Adjustable timer (10-60 minutes)
- 2 year warranty

- 12 modes
- Dual outputs channels
- Rechargeable internal lithium battery
- 20 intensity level settings per channel

Mode description

We suggest that you initially experiment using each of the 12 modes. The mode or combination of modes that gives you the most desirable sensations and comfort is the most appropriate one to use for your current condition.

- Intensity should be set at a level where you will experience some muscular vibration and involuntary muscle movements. It should be set to give a "comfortable" intensity and should not be painful.
- Timing of application should be at least 10 to 20 minutes.
- Frequency of application should be at one to three times per day.

-6-

Directions for use

Three key points for operation:

1. **Find the exact pain point:** The points which the electrode pads are attached to should be the most painful in normal time.
2. **Intensity:** The intensity must be gradually increased, and it's better to reach the highest intensity you can stand.
3. **Application duration:** 2-3 times/day, 40-60 minutes/time, over 10 days for one course.

- Connect the wire to either channel A or B at the top of the unit.
- Connect the pads to the wire by snapping them on.
- Remove the protective film and place the pads on a muscle making sure that both pads are on the skin and not overlapping.
- Turn on the unit via the on/off switch located at the bottom of the unit.
- Select a channel by pressing A/B button at the bottom of the unit.
- When you have turned on the unit, it will automatically start in mode 1.
- To change the modes press the M button.
- Once you have selected a mode, gradually increase the intensity by pressing the + button; to decrease the intensity, press the - button.
- To set the time of use, press the T button.
- Sit back and enjoy the deep soothing sensations!

-7-

*Notes:

The XPD 12 unit is very safe, the output intensity increases only by the "+" key. Even if the intensity is increased to the maximum level, it is within the safe range. When user switches the mode, the intensity will automatically go down to the minimum level.

Electrode guidelines

1. Use only the electrodes supplied by the manufacturer, other electrodes may present a risk of unsuitable electrical characteristics with your stimulator.
2. Do not use the electrodes on different people. Each user must have their own set, otherwise, skin reaction or cross contamination may occur.
3. Always turn power off before removing or repositioning the electrode.
4. Wash skin thoroughly, and then dry it before applying the electrodes.
5. Apply the whole surface of the electrodes firmly to the skin. Do not use electrodes that do not stick properly to the skin or only partially stick to the skin.
6. In case of skin redness under the electrodes after a stimulation session, do not start a new session in the same place if skin redness is still evident.
7. Using electrode pads that are too small or incorrectly applied could result in discomfort or skin burns.

-8-

Electrode recommendations

- The adhesive electrodes that are supplied as a standard accessory with your XPD 12 unit. These are long lasting and reusable pads.
- The recommended electrode pads sizes are bigger rectangular pads for large areas, e.g shoulder, waist, abdomen and leg muscles; large round shape pads for smaller areas, such as forearm muscles, joints etc.

Maintenance and care

- 1 Make sure your skin is free from any dirt, oil, or lotions.
- 2 Before applying the pads, suggest spraying HiDow conductor on your fingers and rub them on both pads. This will help the pads maintain their adhesiveness.
- 3 When you finish using the device, turn it off.
- 4 Take off the pads.
- 5 Place the protective films back on the pads or place the pads on the pad holder.
- 6 Take the wire out from the device and place back in the box until next use.
- 7 Store the XPD 12 unit in a cool, dry place until next use.

- 9 -

Environmental specifications

Operating conditions:

Temperature: 5°C to 40°C (41°F to 104°F)
Humidity (non-condensing) : 30% to 75%
Atmospheric pressure: 700 to 1,060 HPa

Storage and transportation conditions:

Temperature: -10°C to 55°C (14°F to 131°F)
Humidity (non-condensing) : 10% to 90%
Atmospheric pressure: 700 to 1,060 HPa

Technical specifications

Dimensions: 86mm x 52mm x 11mm
Net weight: 52g
Power supply: DC 3.7V

- 10 -

Pulse parameters

XPD 12		TABLE: PULSE parameters		
Rated Supply Voltage (V)	DC 3.7 V			
Load Resistance (Ω)	1000 Ω			
PULSE DURATION (μs)	PULSE repetition frequencies (Hz)	PULSE repetition amplitudes (V)	DC component (V)	Remarks
100	54	50-90	--	Mode 1
100	6	50-90	--	Mode 2
100	1.28	50-90	--	Mode 3
100	54	50-90	--	Mode 4
100	64	50-90	--	Mode 5
100	55	50-90	--	Mode 6
100	54	50-90	--	Mode 7: 7.2
100	<71	50-90	--	Mode 8.1; 8.2; 8.3; 8.4;
100	<60	50-90	--	Mode 9.1; 9.2; 9.3;
100	<48	50-90	--	Mode 10.1; 10.2;
100	8.7	50-90	--	Mode 11
100	<53	50-90	--	Mode 12.1; 12.2; 12.3; 12.4;

- 11 -

Mode description when used as TENS (Pain Relief)

Mode 3: Tapping: simulates a second heart beat. By tapping the muscle you increase the body's natural blood flow and oxygen supply to target area. This also helps flush any inflammation or swelling in target area.

Mode 4: Hot Stone: relaxation mode. This simulates a regular manual rolling massage helping the soft tissue to release and relax.

Mode 6: Shuffle: This is a combination mode. This is a great mode to use when cooling down from an exercise.

Mode 8: Ramping wave: This mode is designed to start out with number 3, tapping. It will get faster and deeper into the belly of the muscle as it continues until it goes to a fully facilitated contraction. Then it will allow the muscle to release. This has been used to increase endurance and build fatigue resistance.

Mode 10: Reflexology: This mode has been designed for the shoes. It stimulates a reflexology session through the socks and/or shoes using several different stimulation patterns.

Mode 11: Auricular Therapy: This mode is also known as "alpha stim" or "micro current". When used properly you should use the earflap attachment. It works 2 ways. By simulating an alpha wave in the brain it has been known to help sleep patterns, behavior, and endorphin release. Also, it will relax the sub occipital muscles and vaso dilate, helping with migraines and headaches.

- 12 -

Mode description when used as EMS (Muscle Stimulation)

Mode 1: Kneading: simulates the thumbs fingers and palms of a real massage therapist. The sensation should feel like a vibration. Focuses on hypertension and stress and helps to break that down returning full range of motion and flexibility.

Mode 2: Acupressure: a much more rapid and focused pulse that focuses on any knot or adhesion in the body. It puts enough pressure on it to break it down and release any acidic (lactic) content within.

Mode 5: Cupping: This will reach in, grab, and stretch the muscle. Help to lengthen and broaden it, the release it as it sets it back. This will help to increase flexibility and mobility of the muscle.

Mode 7: Contract and Hold: This mode is the longest contraction on the machine. It is designed to simulate a flex movement to help build muscle tissue.

Mode 9: Tapping and kneading: a mixture of the fast and slow contractions.

Mode 12: Atlas: It is used for muscle confusion and athletic enhancement.

- 13 -

Indications for use

TENS:

To be used for temporary relief of pain associated with sore and aching muscles in the shoulders, waist, back, neck, upper extremities (arms) and lower extremities (legs) due to strain from exercise or normal household work activities.

Powered Muscle Stimulator (EMS):

It is intended to stimulate healthy muscles including abdomen muscles to improve or facilitate muscle performance.

Great for:

- Stiff neck
- Sore and tight muscles
- Back pain
- Muscular performance
- Aching joints
- Pain relief

Recommended usage points:

Shoulders, neck, shoulder blades, along spine, middle back, lower back, abdomen, Arms: biceps, triceps, forearms, Legs: calves, thighs, Feet: toes, heel, midfoot, joints: shoulder, elbow, wrist, hip, knee, ankle etc.

- 14 -

Unique benefits

Quick acting - Although the rate and degree of response varies for individuals, many users report experiencing positive results after 10 (ten) to 20 (twenty) minutes of use.

No side effects - Since there are no drugs or chemicals involved, HiDow muscle stimulators are safe for long-term use. It is also safe to use while on medications.

Portable - Its pocket travel size enables you to receive treatment anywhere at your convenience.

Cost effective - The device will save you time and money as well as improve your function and enjoyment of life.

Simple to use - High-tech components automate and program most functions. Highly simplified in its functional application.

Be responsible - With the HiDow muscle stimulators, you can take control of your health and well being.

Help others - Everyone suffers from muscle tension in some form. Now you can help others help themselves. Muscle aches and pains may result from working or playing. The problem may be chronic or acute. Everyone, whether they are young or old needs relief from muscle tension and will surely appreciate your gift of health.

Always available - It's always there when you need it: 24 hours a day, every day. You can use it as often as needed.

Extremely effective - The advanced microchip sends out multiple preset electrical, therapeutic, electronic waves through the electrodes. These therapeutic, electronic waves use deep tissue pain relief techniques and various muscle movements. The result is less muscle tension.

Amazingly effective!

- 15 -

Regular application principles

- Find the exact pain point. The points which the electrode pads are attached to should be the most painful point in normal time.
- Intensity - the intensity must be gradually increased and it's better to reach the highest intensity you can stand.
- Application duration: 2-3 times/day, 40-60 minutes/time, over 10 days for one course.
- There are two ways to place the pads, in twin and opposed modes.

Figure 1: Twin mode

Fix the two pads at the top and bottom or both ends on the same side of the human body



Figure 2: Opposed mode

Fix the two pads respectively on the opposite sides of the application position, as shown below

- 16 -

Regular application methods

→ Pain in the trunk

1. Probably pain points in the trunk

According to traditional Chinese medicine (TCM), the most sensible pressure pain point is the key point, the most proper position to be applied. Find the pressure pain point with reference to Fig. 1. Apply one pad onto it. Apply the other pad onto a proper point symmetrical to the chosen pressure pain point or near to it, but the two pads should not be overlapping. If there are several separated pressure pain points, apply them one by one.

2. Recommended positions of pads for application on the neck, shoulder and back (See Fig.3)

Figure 3

Mode 1 or 2 for 30 minutes, Mode 3, 4, 5 or 6 for 30 minutes by turns for A, B, and C. Longer time is needed for those parts of great pain. It is advisable to use fairly large intensity.

3. Recommended positions of pads for application of the waist (See Fig. 4)

Figure 4

Mode 1 or 2 for 30 minutes, Mode 3 or 4 for 20 minutes or Mode 7 or 8 for 30 minutes. It is advisable to do some waist exercise after the application.

- 17 -

→ Pain in the joints and limbs

Pain in the joints, sport injuries, soft tissue injuries, may cause joint aches, which often occur in the joints of shoulder, elbow, knee, wrist, ankle, etc.

1 Position of pads
For joint pain, two pads should be applied to the two symmetrical sides of the painful point, see Fig. 2 (opposed mode).

2 Mode and length of operation time
Mode 1 or 2 for 30 minutes, the Mode 5 or 6 for 30 minutes, 2-3 times a day.

Pain in the muscles and tendons of the limbs, it is often caused by sport injury, overstrain or some other factors. The common signs are muscular pain, swelling and spasm.

- 1 Position of pads
Apply the pads to each end of the painful muscle or apply one pad to the muscle and the other to the tendon. The position can be slightly regulated according to sensation.
- 2 Mode and length of time
Mode 3 or 4 for 20 minutes, then Mode 5 or 6 for 20 minutes. Consult the doctor if the injury has laceration of muscle or tendon.

- 18 -

Electrode placement variations

→ Shoulders

→ Waist and leg

- 19 -

Troubleshooting

Problem	Cause	Solution
The unit is not turning on or I cannot see anything on the screen.	Battery is depleted.	Make sure the unit has been fully charged.
The unit does not seem to be charging.	Charging cable not plugged in fully.	Ensure the USB cable is firmly in the unit on one side and plugged into the AC adapter into the electrical socket on the other side.
Unit does not turn on or the screen is dim.	Battery is low.	Recharge the battery.
Unit is ON, but no sensation is felt from the pads.	Both pads are not firmly on the skin/ There is a loose connection/ The intensity level needs to be increased.	Ensure you are using 2 pads per lead and they are in close proximity to each other, i.e. not on different limbs or too far apart— Make sure both pads are firmly pressed to the skin— Make sure all connections are secure from the unit to the wires and the wires to the pads. Increase the intensity level.

- 20 -

Problem	Cause	Solution
Adhesive electrode pads do not stick to skin even after cleaning and moistening their surface.	Adhesive gel pads need to be replaced.	Replace the set of adhesive gel pads. Pads can be purchased from your local distributor or from www.hidow.com
One pad feels stronger than the other.	This is normal. Different areas will react differently.	Nothing needs to be done. Make sure the pads are moist and are making good contact.
During use the skin feels a painful burning sensation or the stimulation becomes weakened.	Adhesive gel pads are not adhering firmly to the skin. Or the gel pads are too dry.	Put a few drops of water on your fingers and rub them on both pads, which will help the electrodes maintain their adhesiveness. And make sure the pads are pressed firmly to the skin during application.
Did not provide adequate relief.	Not using it long enough or frequently enough.	Use your XPD 12 20-60 minutes at a time, 3-6 times per day. Recommended: lying down
Only 1 port on the unit works.	Damaged port or poorly conducting pads or lead wire.	Isolate the problem by trying another lead and set of pads that you know that work in the troublesome port.

- 21 -

Contraindications

- DO NOT use this device if you have any of the following medical conditions:
 1. An implanted cardiac pacemaker, defibrillator or other implanted metallic or electronic device
 2. Undiagnosed pain syndromes
 3. Have been diagnosed with cancer
 4. Are pregnant
 5. Have suffered acute trauma or surgical procedure in the past six months
 6. Have cardiac problems or cardiac disease
 7. Have epilepsy
 8. Have painful and/or atrophied muscles
 9. Have abdominal or inguinal hernia
 10. Have limited range of motion in skeletal joints
 11. Have blood circulatory problems
- The XPD 12 unit is intended for use by ADULT and ONLY on healthy muscles. NEVER use this device on muscles that are atrophied, painful, suffer spasms, on a limb with painful or otherwise afflicted joints.
- The XPD 12 unit is not intended for the application of any medical condition or disease nor is it intended for physiotherapy or muscle rehabilitation. It is contraindicated for use on any muscle that is injured or diseased.

- 22 -

- Do not attempt to use XPD 12 unit for muscle reeducation, to prevent muscle atrophy or spasms, improving range of motion, blood flow deficiencies/venous thrombosis.

Warnings

- 1 Stimulation should not be applied over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
- 2 Stimulation should not be applied transthoracically in that the introduction of electrical current into the heart may cause cardiac arrhythmias.
- 3 Stimulation should not be applied transcranially.
- 4 Apply electrode ONLY to normal, intact, clean skin. Do not apply electrodes over open wounds or over swollen, infected or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, varicose veins, etc.
- 5 DO NOT APPLY STIMULATION:
 - Over frontal area of the neck (near site of carotid sinus nerves)
 - Over the neck or mouth. Severe spasms of the laryngeal and pharyngeal muscles may occur with contractions strong enough to and/or cause difficulty in breathing. Stimulation over the neck could also have adverse effects on the heart rhythm or blood pressure.
 - Transcranially.
 - Over the swollen, infected or inflamed areas of skin eruptions (e.g. phlebitis, thrombophlebitis, varicose veins, etc).

- 23 -

- Across the chest. Consult your physician before using this device because it may be possible to cause lethal rhythm disturbances to the heart in susceptible individuals.

- Over, or in proximity to cancerous lesion.

- If you are epileptic

- After experiencing acute trauma or fracture

- Following recent surgery

- If you have a hernia (abdominal or inguinal)

- To the frontal, laryngeal and temporal region of the neck.

- 6 Never use the XPD 12 unit while driving, operating machinery or during activities in which involuntary muscle contractions may endanger the user or others.

- 7 The effects of stimulation of the brain are unknown. Therefore, do not apply stimulation across the head and do not place electrodes opposite sides of the head.

- 8 Do not use the XPD 12 unit in the bath or shower.

- 9 Persons with suspected heart problems or epilepsy should obtain opposite medical advice.

- 10 Never use the XPD 12 unit while sleeping.

- 11 Never immerse the XPD 12 unit in any liquid.

- 24 -

Precautions

- This XPD 12 unit should be used only with the leads, electrodes and accessories provided by the manufacturer.
- This XPD 12 unit should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.

Adverse reactions

- You may experience skin irritation and burns beneath the stimulation electrodes applied to your skin;
- You may experience headache and other painful sensations during or following the application of electrical stimulation near your eyes and to your head and face; and
- You should stop using the XPD 12 unit and should consult with your physician if you experience adverse reactions from the unit.

Location of stimulation

- Stimulation near metal:** Do not use the stimulator when close to metal. This increases the potential for discomfort and burns. Remove jewelry, piercings, belt clips or any other removable metal objects in the area of stimulation.
- Stimulation of eyes:** Do not use the stimulator directly to the eyes.
- Stimulation of the neck or mouth:** Do not use the stimulator on the neck (particularly the carotid sinus or mouth). This could cause severe muscular spasms causing airways to close, breathing difficulties, or side effects to heart rhythm or blood pressure.

- 25 -

Stimulation of healthy skin: Only use the stimulator on healthy skin. Stimulation can worsen disorders caused by lesions or skin diseases.

Stimulation across the torso: Do not use the stimulator across the torso (one electrode placed on the torso and one electrode placed on the back). The electrical current passing through the torso can cause potentially fatal heart rhythm disruptions.

Cranial stimulation: The effects of stimulation on the brain are not known. For this reason stimulation should not occur to the head, and electrodes must not be placed either side of the head.

Conformity to safety standards

Statement of EMC
IEC 60601-1-2:2014

Conformity to MDD Requirements

IEC 60601-1:2005+AMD1:2012/ANSI/AAMI ES60601-1:2005+A1:2012, IEC 60601-1-11:2015, IEC 60601-2-10:2012

- 26 -

Warranty information

- Warranty covers defects in materials and workmanship under normal use for a period of 2 years from the date of retail purchase by the original end-user purchaser ("Warranty Period"). If a defect arises and a valid claim is received within the Warranty Period Hi-Dow will either repair the defective unit or exchange the product with a product that is new or which has been manufactured from new or serviceable used parts.
- To register your product warranty please visit www.hidow.com, click on the warranties menu and select the warranty that refers to your product: standard 2-year or lifetime VIP.

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