



dpl[®]

Oral Care

Light Therapy System

INCLUDED

- Oral Care Control Device
- Interchangeable Mouth Tray
- Power Supply
- Power Cable
- Storage Case
- Silicone Cleaning Cover
(Do Not Discard)



Interchangeable Mouth Tray



Power Supply

Power Cable

Oral Care Control Device

WELCOME

Welcome to LED Technologies, Inc and our dpl® Oral Care line. We are thrilled to announce this product as an in-home device to be used at your convenience.

The dpl® Oral Care Light Therapy System uses light emitting diodes (LEDs) in the infrared and red spectrum, which are medically proven to **increase circulation, enhance gum health, and relieve pain**. Non-UV blue lights **increase the effectiveness of whitening gels by removing tannic stains from coffee, tea and wine, and leaving whiter teeth**.

Please read this manual before operating the Oral Care System.

CHARGING INSTRUCTIONS

- **Only use provided Power Supply. Attempts to use a third-party power supply could damage device and will void the warranty.**
- **Charging:**
 - **Before first use, charge Oral Care Control Device.**
 - **Charging typically takes 4-6 hours but may take longer.**
 - Separate the Oral Care Control Device from the Interchangeable Mouth Tray.
 - Plug Power Cable into Power Supply and Oral Care Control Device.
 - Plug Power Supply into power source.
 - While charging, red light flashes on Oral Care Control Device on/off button.

- When fully charged, red light is solid red on the Oral Care Control Device on/off button.
- **Battery Indicator:**
 - The Battery Indicator Light displays on the On/Off Button.
 - Check Oral Care Control Device On/Off button:
 - If the button is **green**, device is charged and ready for use.
 - If the button is **red**, device needs to be charged.
 - LED lights on the Mouth Guard Tray light up when On/Off button on the Oral Care Control Device is pressed.
 - The Oral Care Gum Device will automatically shut off after 10 minutes.
 - This time may vary based on the battery charge level of the Oral Care Control Device.
 - If automatic shut off occurs quicker than 10 minutes, charge Oral Care Control Device.

OPERATING YOUR ORAL CARE SYSTEM

- To turn device On, press and hold the On/Off button for 2 seconds or until LED lights on the Interchangeable Mouth Tray light up.
- Place device in mouth.

TREATMENT

FOR TEETH WHITENING:

- Treatment time is 20 minutes. Device will automatically shut off after 10 minutes. (Requires two 10 minute sessions.)

- Use with teeth whitening gel of your choice.
- For first week, depending on your whitening gel strength and teeth sensitivity, use once or twice daily (treatments can be back-to-back).
- After first week, use 2-3 times per week or as needed.
- Helpful hint for teeth whitening: for maximum results, put gel directly on teeth instead of on tray.

FOR GUM HEALTH:

- Use once daily for a minimum of 10 minutes.

Troubleshooting:

1. System is not turning on:
 - Is the device fully charged? (Power indicator lights will be on and solid green.)
 - If device still does not turn on, call Customer Service at 1 (800) 337-9565.
2. System is not charging:
 - Is USB charging cord firmly inserted into device? (Does power indicator light turn on?)
 - Try different USB power adaptor to ensure it is working properly.
 - Try different wall outlet to determine if there is power to the outlet.
 - If device still does not charge, call Customer Service at 1(800) 337-9565.

HOW TO CLEAN DEVICE

- While washing, pull apart device and install waterproof Silicone Cleaning Cover over Interchangeable Mouth Tray's USB connection.
- Alternatively, you may carefully cover the USB plug with your thumb to avoid contact of USB plug with water.

CONNECT WITH US

When you register your **dpl**[®] Oral Care Light Therapy System at revivelighttherapy.com/warranty, you will be entered into our system to be eligible for exclusive offers, VIP discounts, new product previews, and the opportunity to earn free products through new product trials.

DISCLAIMER

This product is not intended to diagnose, cure or prevent disease. LED Technologies, Inc makes no claims, representations or warranties regarding the ability of this product to cure any physical, oral or mental condition. A qualified health professional should be consulted with regard to any condition requiring medical attention.

WARRANTY

LED Technologies, Inc warrants that this product will be free from defects in materials and workmanship for 6 months from the date of purchase by the original owner. This warranty covers the original purchaser only and is not assignable or transferable to another party. Register your warranty online at revivelighttherapy.com/warranty



Our products are continually improving. The actual product may differ from the product shown.