

CITY COLLECTION

Battery Heated Gloves User's Guide

Step 1: Fully Charge the Batteries Plug the adaptor into the wall and connect the batteries using the splitter. The **Red** light on the adaptor indicates that the batteries are charging. The **Green** light indicates that the batteries are charged.

Step 2: Connect the Battery On the back of the gloves, locate the zippered pocket and find the plug inside. Connect the plug to the batteries. Notice that the **+** button on the top of the glove will flash yellow, then turn off. This indicates that the batteries are properly connected. Next, place the batteries inside of the zippered pocket on either side of the plug.



Step 3: Power On Press and hold the **+** button on the top of the glove for 2 seconds until the button lights up in red.

Step 4: Adjust the Temperature Click the + button on top of the gloves to cycle through the 3 temperature settings. (*Heating time may vary based on condition of the battery)

Red = High Heat (2Hrs) / Yellow = Medium Heat (4.5 Hrs) / Green = Low Heat (8 Hrs)

Step 5: Power Off Press and hold the **+** button on the top of the glove for 2 seconds until the light on the button turns off.

Storage Tips Store in a cool dry place. Fully charge the battery before each use. Fully charge the battery prior to storing and at least once every 6 months of storage thereafter.

Wash and Care Remove the battery from the garment prior to washing. Hand wash in cold water only. Hang dry. Do not iron. Do not use bleach or cleaning solvents. Do not twist or wring. Do not turn power on when garment is folded.

Warning Avoid direct contact with all liquids. Never tamper with connection or use other power source.

Helpful Reminders

Do not apply the garments to sensitive skin areas, acute injuries, and open wounds. This garment is not recommended for use by children, pregnant women, and immobile individuals with disabilities. In the unlikely occurrence of skin discomfort and irritation, immediately discontinue the use of the product and consult your physician and health care professional.
