



Posey® Incontinence Bed Pad

BED CANOPY SYSTEM

INCONTINENCE BED PAD

The Posey Incontinence Bed Pad is designed for use inside the Posey Bed Canopy System, or whenever an incontinent bed pad is needed. The 8005A is designed with pockets to fit the 8007 Support Surface. The multi-layer construction ensures that the Posey Incontinence Bed Pad will absorb moisture, without penetrating through the backing and onto the patient's bedding. Provides superior absorption and wicking properties, to protect both the patient's skin and bedding from urine. Machine washable. 35"x 82". One per package; Case qty. of two.

REF 8005A Incontinence Bed Pad for Support Surface

APPLICATION INSTRUCTIONS:

The Posey Incontinence Bed Pad should be positioned under the patient with the blue side facing down. To secure pad in place, place each pocket on the opposite ends of the pad onto each end of the 8007 Support Surface. Once the pad is in place, check to ensure that the patient cannot become trapped under the pad.



WASHING INSTRUCTIONS:

Wash lukewarm, bleach as directed, dry on low.



J.T. Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA
Phone: 800-447-6739 • Fax: 800-767-3933 • Web: www.posey.com



MDSS
Burckhardtstr. 1,
D-30163, Hannover, Germany
© 2006 J.T. Posey Company. All rights reserved.



M6190 081106



Posey® Incontinence Bed Pad

BED CANOPY SYSTEM

INCONTINENCE BED PAD

The Posey Incontinence Bed Pad is designed for use inside the Posey Bed Canopy System, or whenever an incontinent bed pad is needed. The 8005A is designed with pockets to fit the 8007 Support Surface. The multi-layer construction ensures that the Posey Incontinence Bed Pad will absorb moisture, without penetrating through the backing and onto the patient's bedding. Provides superior absorption and wicking properties, to protect both the patient's skin and bedding from urine. Machine washable. 35"x 82". One per package; Case qty. of two.

REF 8005A Incontinence Bed Pad for Support Surface

APPLICATION INSTRUCTIONS:

The Posey Incontinence Bed Pad should be positioned under the patient with the blue side facing down. To secure pad in place, place each pocket on the opposite ends of the pad onto each end of the 8007 Support Surface. Once the pad is in place, check to ensure that the patient cannot become trapped under the pad.



WASHING INSTRUCTIONS:

Wash lukewarm, bleach as directed, dry on low.



J.T. Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA
Phone: 800-447-6739 • Fax: 800-767-3933 • Web: www.posey.com



MDSS
Burckhardtstr. 1,
D-30163, Hannover, Germany
© 2006 J.T. Posey Company. All rights reserved.



M6190 081106