



Indication:

Patients at high risk for pressure sores and/or requiring posture support and/or inclined to slide forward.

Patient Benefits:

The Posey Molding Foam Cushion is constructed with 1½ thick viscoelastic foam (“Memory foam”), laminated to a firm urethane foam base. This combination provides the benefits of stability and even pressure distribution without constricting blood circulation, lessening the risk of decubitus ulcer development.

- Unique “rounded edge” design minimizes pressure under the knees and thighs
- Convex shape reduces the “hammocking” effect of sling wheelchair seats; provides better stability and pressure equalization
- Covered in Posey’s exclusive Sure-Grip™ vinyl
 - Waterproof
 - Anti-microbial to protect the fabric
 - Four-way stretch vinyl promotes more effective immersion into the molding foam
 - Sure-Grip™ surface prevents sliding and helps hold the cushion in place
- Zippered outer cover removes quickly and easily for laundering
- Meet CA #117 flame retardancy standards
- Optional “PoseyAire” fabric cover with ‘safety straps’ for secure chair attachment and superior moisture wicking capabilities
- Posey cushions are provided in convenient box packaging with “punch in” handle for easier storage and product identification



Unique memory foam provides superior immersion compared to standard urethane foam.

Laundering Instructions:

Vinyl covers wipe clean with disinfectant, or can be removed for machine washing.



WIPE CLEAN

OR



WASH HOT
25 MIN.



BLEACH AS DIRECTED
ON CONTAINER



DRY ON
LOW

Warranty:

The Posey Molding Foam cushion is warranted to perform as indicated. The Posey Company will, at their option, repair or replace, except for vandalism or other abuse, any cushion for a period of one (1) year from date of purchase.

CUSHION	DIMENSIONS W"xD"xHb-Hf"	BOTTOM DESIGN	WEIGHT	COVER MATERIAL	COVER COLOR	POSEYAIRE COVER
7024C	18 x 16 x 3	Convex	3 lbs.	Sure-Grip™	Dk. Blue	7422