

## **Relieve Tight Muscles and Pinched Nerves**

The **ChiSoft** cervical neck traction device is indicated for tight muscles, pinched nerves, neck pain, tension headaches, osteoarthritis, swollen discs, straightening of the neck and for degenerative discs or joints. The user should consult a medical professional if preexisting conditions exist.

**ChiSoft** Neck Traction is comfortable and can be used anywhere. The **ChiSoft** Cervical Neck Traction device is an air-inflated neck-stretching device that treats most neck problems. This Air Neck Traction device is made up of lightweight soft rubber and covered with soft textile at the top and bottom, remaining portable and convenient for users on the go.

### **Who can benefit from **ChiSoft** Neck Traction?**

This portable Air Neck Traction device has several built-in air chambers and provides an effective, gentle, even traction particularly suitable for office workers and people with whiplash. For the Air Traction Devices, use the inflation pump to slowly raise the pressure to a comfortable stretch. The traction allows the tight neck muscles to relax and relieve the pressure in the joints and nerves.

### **Cervical Neck Traction can be used anywhere**

This portable Neck Traction Comforter has several built-in air chambers and provides an effective, gentle, even traction particularly suitable for office workers and people with whiplash. For the Air Traction Devices, use the inflation pump to slowly raise the pressure to a comfortable stretch. The traction allows the tight neck muscles to relax and relieve the pressure in the joints and nerves.

## How to use Cervical Traction Device?

Simply place the Neck Traction device around your neck to bring relief on the spot. This cervical Neck Traction device can be taken anywhere to massage your neck. No assembly or batteries are required to operate the unit. The simplicity of the design makes it the ideal tool for home or office use. No longer do you have to wait until you get home to seek relief, or lie on the floor to use another product.


### Specifications - ChiSoft Neck Traction Deluxe

#### Safety & Precautions

- Inner Diameter: **12 cm**
- Outer Diameter: **27cm**
- Neck Size: **11" to 19"**
- Weight: **8.8 oz (250 gram)**
- Inflatable Layers: **3**
- Adjustable Height: **YES**
- Material: **Velour**
- Max Rising High: **21cm / 8.3inch**



- Do not sleep with the Air Neck Traction
- Do not use as neck brace
- Do not use with acute neck strains or sprains.
- Adjustable Size
- Adjust the cushion vertically by blocking the upper layer to be inflated



### **Special Feature ChiSoft Neck Traction PROII?**

More comfortable stretch at the front of the neck.

