



**INDICATION:**

Posey's Wedge Foam Pommel Cushions are designed to keep hips and knees in alignment, while minimizing patient forward sliding. Can be used to help prevent hip dislocation following surgery and allow patients to exercise their knees and lower legs. Ideal for use in standard wheelchairs or geri-chairs.

**PATIENT BENEFITS:**

- Foam inner construction is lightweight and provides exceptional comfort and support.
- Both Sure-Grip and standard vinyl material used on all Posey cushions are warranted to be waterproof and anti-microbial to protect the fabric.
- Built-in pommel assists hip and knee alignment.
- Zippered outer cover can be quickly removed for cleaning.
- Posey's exclusive Sure-Grip vinyl provides a stable, comfortable seating surface.
- Non-skid bottom to help prevent sliding on seating surfaces.
- Convex bottom design can be used in "sling wheelchairs" without the need for an additional foundation.
- Meets CA #117 flame retardancy standards.


**CLEANING INSTRUCTIONS:**

Vinyl covers wipe clean with mild detergent and/or disinfectant, or can be removed for machine washing.



WIPE CLEAN

OR


 WASH HOT  
25 MIN.

 BLEACH AS DIRECTED  
ON CONTAINER

 DRY ON  
LOW

**WARRANTY:**

All Posey Wedge Foam Cushions are warranted to perform as indicated. The Posey Company will, at their option, repair or replace, except for vandalism or other abuse, any cushion for a period of one (1) year from date of purchase.

CAT. #	BOTTOM	DIMENSIONS	WEIGHT
<b>6820F</b>	Flat	18" x 16" x 3 1/2" - 1 1/2"	2 lbs.
<b>6820C</b>	Convex	18" x 16" x 3 1/2" - 1 1/2"	2 lbs.
<b>6830C</b>	Convex	20" x 16" x 3 1/2" - 1 1/2"	2.5 lbs.
<b>6830F</b>	Flat	20" x 16" x 3 1/2" - 1 1/2"	2.5 lbs.

Pommel measures 4"W x 5"D x 4"H