

Purpose


Skil-Care's Econo-Gel Pads are designed to relieve pressure and increase comfort for wheelchair users. High-density foam layers on top and bottom provide stability and cushioning. The gel-filled center inserts add pressure-relieving flotation.

Important Notice

Econo-Gel Pads have safety ties that must be secured to the wheelchair to prevent the cushion from slipping off the wheelchair seat.

Easy Care

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser or mild disinfectant. The polyester cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.



Econo-Gel Pads - Standard Series

Re-Order Number	Cover	Size
751150	Vinyl w/Polyester	16" x 16" x 2"
751151	Vinyl w/Polyester	18" x 16" x 2"
751148	Vinyl w/Polyester	18" x 18" x 2"
751149	Vinyl w/Polyester	20" x 18" x 2"
751153	Vinyl w/Polyester	20" x 16" x 2"
751154	Vinyl Only	16" x 16" x 2"
751155	Vinyl Only	18" x 16" x 2"
751157	Vinyl Only	20" x 16" x 2"
751156	Vinyl Only	18" x 18" x 2"




Econo-Gel Pads - Bimini Series

Re-Order Number	Cover	Size
751180	Bimini Blue	16" x 16" x 2"
751181	Bimini Blue	18" x 16" x 2"
751184	Bimini Pink	16" x 16" x 2"
751185	Bimini Pink	18" x 16" x 2"

Econo-Gel Pads - Thin-Line Series

Re-Order Number	Cover	Size
751170	Vinyl w/Polyester	16" x 16" x 1.5"
751171	Vinyl w/Polyester	18" x 16" x 1.5"
751174	Vinyl w/Polyester	20" x 16" x 1.5"
751178	Vinyl w/Polyester	18" x 18" x 1.5"
751179	Vinyl w/Polyester	20" x 18" x 1.5"
751172	Vinyl Only	16" x 16" x 1.5"
751173	Vinyl Only	18" x 16" x 1.5"
751175	Vinyl Only	20" x 16" x 1.5"
751176	Vinyl Only	18" x 18" x 1.5"
751177	Vinyl Only	20" x 18" x 1.5"



Turn product insert over for more of Skil-Care's Pressure Control Cushions, Wheelchair Accessories and Bed Positioning products.