

Re-Order #	Description	Size	Units
754750	Sling-Seat Cushion w/LSII Cover	16"W x 16"D x 3"H	EA
754755	Sling-Seat Cushion w/LSII Cover	18"W x 16"D x 3"H	EA
754757	Sling Seat Cushion w/LSII Cover	18"W x 18"D x 3"H	EA
754810	Sling Seat Cushion w/LSII Cover	20"W x 16"D x 3"H	EA
754812	Sling Seat Cushion w/LSII Cover	20"W x 18"D x 3"H	EA
754814	Sling Seat Cushion w/LSII Cover	22"W x 16"D x 4"H	EA
754816	Sling Seat Cushion w/LSII Cover	22"W x 18"D x 4"H	EA
754818	Sling Seat Cushion w/LSII Cover	24"W x 18"D x 4"H	EA
754821	Sling Seat Cushion w/LSII Cover	24"W x 20"D x 4"H	EA
754823	Sling Seat Cushion w/LSII Cover	26"W x 18"D x 4"H	EA
754825	Sling Seat Cushion w/LSII Cover	26"W x 20"D x 4"H	EA
754827	Sling Seat Cushion w/LSII Cover	28"W x 18"D x 4"H	EA
754829	Sling Seat Cushion w/LSII Cover	28"W x 20"D x 4"H	EA
754833	Solid Seat Cushion w/LSII Cover	20"W x 16"D x 3"H	EA
754835	Solid Seat Cushion w/LSII Cover	20"W x 18"D x 3"H	EA
754837	Solid Seat Cushion w/LSII Cover	22"W x 16"D x 4"H	EA
754839	Solid Seat Cushion w/LSII Cover	22"W x 18"D x 4"H	EA
754841	Solid Seat Cushion w/LSII Cover	24"W x 18"D x 4"H	EA
754843	Solid Seat Cushion w/LSII Cover	24"W x 20"D x 4"H	EA
754845	Solid Seat Cushion w/LSII Cover	26"W x 18"D x 4"H	EA
754847	Solid Seat Cushion w/LSII Cover	26"W x 20"D x 4"H	EA
754849	Solid Seat Cushion w/LSII Cover	28"W x 18"D x 4"H	EA
754851	Solid Seat Cushion w/LSII Cover	28"W x 20"D x 4"H	EA



Custom sizes are available. Please call Skil-Care at 800-431-2972.

Purpose

Skil-Care's Tri-Foam Cushions provide maximum comfort and pressure reduction for bariatric wheelchair users. The three-part design incorporates a top layer of temperature sensitive visco-elastic foam that molds to the user to eliminate pressure "hot spots." The center layer of high-density foam offers comfort cushioning, and protection against bottoming-out. The firm base distributes the user's weight for efficient pressure reduction.

Low-Shear II cover fabric

This liquid-resistant, double-stretch, breathable polyurethane fabric with laminated polyester backing reduces skin-damaging shear and friction.

Important Safety Notice

Be certain to secure the cushion's safety strap to the wheelchair frame to prevent the cushion from sliding off the seat.

Easy Care

The Low-Shear II cover wipes clean with a damp cloth. Mild spray-type cleansers and disinfectants may be used. Do not launder the cover.