







13604*



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Pull Cord Alarm

13602 1/ea

13615 Replacement Cord with Clip, 1/ea

(For use with 13602)

FEATURES:

- Economical way to notify staff that a patient has moved.
- Alarm easily secures to a bed or chair with a self contained clip.
- Activation cord is adjustable from 18" to 36".
- When user gets up, the pin connector pulls from the alarm unit activating the signal.
- Deactivate with On/Off switch.
- Two volume settings: 85 dB and 100 dB.
- Requires 1 9V Battery (Included).
- 1 Year Limited Warranty.

New & Improved Tamper Proof Magnetic Pull Cord Alarm

(Comes with Screwdriver to Remove Battery Cover)

13603 1/ea

13617 Replacement Cord, Magnet and Clip, 1/ea

(For use with 13603)

13618 Replacement Plastic Alarm Caddy, 1/ea

FEATURES:

- When magnetic pull switch is removed, alarm is activated.
- Pull switch can be activated from any direction.
- Convenient On/Off switch for easy activation-deactivation.
- Easily mounted on wheelchair or bed.
- Cord Adjusts from 28" to 58".
- Alligator clip allows for attachment to patient's clothes.
- 85-95 dB.
- One 9 volt battery included.
- 1 Year Limited Warranty.

Tamper Proof Seat Belt Alarm

13604* 1/ea

FEATURES:

- Alarm sounds when belt is unbuckled.
- Can be set for 85 dB or 100 dB.
- Alarm has convenient On/Off switch.
- LED flashes with the alarm to assist in locating the patient.
- Includes a bracket to mount the alarm box to the wheelchair.
- One 9 volt battery included.
- 1 Year Limited Warranty.

Alarm Door Banner

13098 1/bx

FEATURES:

- High visibility combined with a magnetically activated alarm create an effective deterrent to wandering.
- Arresting yellow banner and brilliant, red "STOP" sign.
- Alarm features an off/soft/loud switch and LED low-battery indicator.
- Attaches to doorjamb with hook-and-loop.
- Fits doors from 30" to 42" (W).
- Requires two AAA batteries, sold separately.
- 6 Month Warranty.