





Designed to treat common ear conditions, such as :

- Otitis Media with Effusion
- Aerotitis/Barotitis (caused by rapid elevation changes)
- Eustachian Tube Dysfunction

Based on the Politzer Maneuver, the EarPopper® delivers a safe, constant, regulated flow of air into the nasal cavity. During the moment of swallowing, the air is diverted up the eustachian tube, clearing and ventilating the middle ear. The EarPopper® relieves negative ear pressure and allows accumulated fluids to drain. Non-sterile. Reusable.



EP-2100 EarPopper®, Home 1 EP-1010 Disposable Tips 50

HOME



**PROFESSIONAL / OFFICE VERSION** 

Prod	duct No.	Description	Units/Box
EP-	3000	EarPopper®, PRO	1
EP-	1000	Disposable Tips	50



Units/Box