

molift[®]
a part of Etac



Molift EvoSling

Comfort, functionality
and safety

etac[®]

Molift EvoSling



Molift has been manufacturing high quality hoists and slings for generations. Quality, comfort and safety are key considerations during product development. We offer a wide range of sling models and sizes ensuring we have a solution, whatever the users requirements.

What are the user's requirements?

For people who need help with mobility, it's important that the sling is both comfortable and safe. It must facilitate transfers and meet the user's requirements.

Comfort, function and safety

Our slings have been developed in conjunction with carers and users, in order to meet their requirements in terms of functionality, comfort and safety. The slings must be comfortable and secure, while at the same time being easy to work with for both the user and carer.

Our standard models are designed for use in the usual types of transfers and are suitable for most users. They are available in three different back heights: HighBack, MediumBack and LowBack. We know that users may have a wide variety of requirements, so we've designed additional models with specific features to meet their needs, including the EvoSling Hygiene and EvoSling Ampu.

The design of our slings ensures that they distribute weight and pressure evenly for maximum safety and comfort.

Individual testing

We recommend individual testing, and offer a wide range of different models and sizes. To assist therapists with this we have a Sling Assessment Guide.



Colour coding

Colour coding of hoisting loops helps you find the right fit.

Padding

Padding in back and leg support distributes the weight over a larger area, making the sling more comfortable for the user.

Belt

EvoSling Hygiene has a belt for extra safety as well as anti-slip protection on the inside.

Guide straps

Guide straps assist with positioning during the lifting process.

EvoSling MediumBack Padded.

Flexibility

Multiple hoisting loops on both head and leg supports allow the distance to the suspension and the user's sitting angle to be adjusted.

Three different back heights

HighBack, MediumBack and LowBack to suit different needs.

Labelling

Clear, practical information on a highly visible label. Washing instructions and technical information are found in a protective pocket at the back.

Application pocket

Pockets in the back and on the leg supports make application easier.

Informative, clear labelling

The large label contains all the information you will need for day-to-day use. Washing instructions, serial number and other technical information can be found behind the label in a pocket which protects against wear and tear during washing.


Periodic inspection

The sling has a label on which you can record when a periodic inspection has taken place. It can be found in the label pocket. By recording the date on which the sling is first used on the label, you can plan when the first inspection should take place.


Materials

The polyester slings are hardwearing and durable. The material is flexible and easy to maintain and has a dirt-repellent surface.

The polyester net slings are equally hardwearing, durable and flexible. The material also offers good ventilation.



Polyester net



Polyester

Periodic inspection performed

| MONTH | YEAR | | | | | | | |
|-------|------|------|------|------|------|------|------|------|
| | 2014 | 2015 | 2016 | 2017 | 2018 | 2019 | 2020 | 2021 |
| 1 | | | | | | | | |
| 2 | | | | | | | | |
| 3 | | | | | | | | |
| 4 | | | | | | | | |
| 5 | | | | | | | | |
| 6 | | | | | | | | |
| 7 | | | | | | | | |
| 8 | | | | | | | | |
| 9 | | | | | | | | |
| 10 | | | | | | | | |
| 11 | | | | | | | | |
| 12 | | | | | | | | |

User name / Navn

.....

.....

.....

Detailed information about EvoSling

Recommended user weight

| | |
|-----|------------|
| XS | 17–25 kg |
| S | 25–50 kg |
| M | 45–95 kg |
| L | 90–160 kg |
| XL | 160–240 kg |
| XXL | 230–300 kg |

Materials

Polyester
PE-cellfoam (padded models)

Washing instructions



Washing at higher temperatures wears out the material faster.

Colour coding



CONTENTS

molift[®]
a part of Etac

| | |
|-------------------------------------|----|
| Molift EvoSling | 2 |
| Molift EvoSling HighBack | 7 |
| Molift EvoSling MediumBack | 8 |
| Molift EvoSling LowBack | 9 |
| Molift EvoSling Hygiene | 10 |
| Molift EvoSling Amпу MediumBack | 11 |
| Molift EvoSling Comfort MediumBack | 13 |
| Accessories for the Molift EvoSling | 14 |

www.etac.com



molift
AWP 205

Molift EvoSling HighBack

FOR USERS WHO NEED A LOT OF SUPPORT

Comfortable support

EvoSling HighBack supports the entire back and has a comfortable head support. The sling is suitable for most users, including those with poor core and head stability.

EvoSling HighBack is suitable for lifting to and from a lying position, as well as from the floor. The sling provides a comfortable reclined sitting position.

The pocket in the lower back and on the leg supports facilitates application of the sling.

Two options

Molift EvoSling HighBack comes with two options: polyester with padding in the back and leg supports and a soft polyester net with padded head and leg supports.




Easier positioning

The polyester net sling has a red seam in the middle of the back to help ensure that it is positioned symmetrically.




Molift EvoSling HighBack



| Size | ① Padded Item no.: | ② Net Padded Item no.: |
|------|-----------------------|---------------------------|
| XS | 1722810 | 1722710 |
| S | 1722820 | 1722720 |
| M | 1722830 | 1722730 |
| L | 1722840 | 1722740 |
| XL | 1722850 | 1722750 |
| XXL | 1722860 | 1722760 |

| Size | Ⓐ Back height | Ⓑ Width |
|------|---------------|---------|
| XS | 70 cm | 30 cm |
| S | 80 cm | 35 cm |
| M | 93 cm | 40 cm |
| L | 100 cm | 47 cm |
| XL | 105 cm | 54 cm |
| XXL | 110 cm | 60 cm |



Molift EvoSling MediumBack

FOR USERS WHO HAVE HEAD CONTROL BUT NEED FULL BACK SUPPORT

Full back support

EvoSling MediumBack is a versatile model which is suitable for the most common lifting situations. The sling provides a comfortable upright sitting position with full back support.


Users place their arms inside the sling. A pocket in the lower back and on the leg supports facilitates application of the sling. The polyester net sling has a red seam in the middle of the back to help ensure that it is positioned symmetrically.

Two options

Molift EvoSling MediumBack comes with two options: polyester with padding in the back and leg supports and a soft polyester net with padded leg supports.




Molift EvoSling MediumBack



| Size | ① Padded Item no.: | ② Net Padded Item no.: |
|------|-----------------------|---------------------------|
| XS | 1723010 | 1723110 |
| S | 1723020 | 1723120 |
| M | 1723030 | 1723130 |
| L | 1723040 | 1723140 |
| XL | 1723050 | 1723150 |
| XXL | 1723060 | 1723160 |

| Size | Ⓐ Back height | Ⓑ Width |
|------|---------------|---------|
| XS | 54 cm | 30 cm |
| S | 56 cm | 37 cm |
| M | 60 cm | 49 cm |
| L | 67 cm | 57 cm |
| XL | 70 cm | 72 cm |
| XXL | 73 cm | 81 cm |



Molift EvoSling LowBack

FOR ACTIVE USERS WITH GOOD CORE STABILITY


Shoulder-height support

EvoSling LowBack is a versatile model which provides shoulder-height support. The sling provides a comfortable sitting position and users can place their arms inside or outside the sling.


A pocket in the lower back and pockets on the leg supports facilitate application of the sling. The polyester net sling has a red seam in the middle of the back to help ensure that it is positioned symmetrically.

Three options

Molift EvoSling LowBack comes with three options: polyester with padding in the back and leg supports and soft polyester net with or without padded leg supports.




Molift EvoSling LowBack



| Size | ① Padded Item no.: | ② Net Item no.: | ③ Net Padded Item no.: |
|------|-----------------------|--------------------|---------------------------|
| XS | 1722310 | 1722410 | 1722510 |
| S | 1722320 | 1722420 | 1722520 |
| M | 1722330 | 1722430 | 1722530 |
| L | 1722340 | 1722440 | 1722540 |
| XL | 1722350 | 1722450 | 1722550 |
| XXL | 1722360 | 1722460 | 1722560 |

| Size | Ⓐ Back height | Ⓑ Width |
|------|---------------|---------|
| XS | 41 cm | 40 cm |
| S | 43 cm | 43 cm |
| M | 44 cm | 48 cm |
| L | 47 cm | 55 cm |
| XL | 51 cm | 62 cm |
| XXL | 55 cm | 70 cm |



Molift EvoSling Hygiene

FOR USERS WITH GOOD CORE STABILITY AND MUSCLE TONE

Easy to position


EvoSling Hygiene has a low back and a large leg support opening to facilitate dressing and undressing when using the toilet.

Application pockets



Molift EvoSling Hygiene is made of polyester and has a padded back and padded leg supports. The application pockets on the leg supports makes positioning them easier.

Extra security

EvoSling Hygiene has a belt for added security during lifting.



Molift EvoSling Hygiene



Size

XS
S
M
L
XL
XXL

Padded Item no.:

1722210
1722220
1722230
1722240
1722250
1722260

Size


XS
S
M
L
XL
XXL

(A) Back height

29 cm
33 cm
33 cm
33 cm
34 cm
35 cm

(B) Width of belt

35 cm
46 cm
52 cm
70 cm
73 cm
76 cm



Molift EvoSling Ampu MediumBack

FOR USERS WITH AN AMPUTATED LEG OR WITH POOR MUSCLE TONE

Secure and comfortable

EvoSling Ampu MediumBack has a narrow leg support opening and a shoulder-height back. The sling's leg supports are wrapped under the legs. The leg support straps are passed through the strap guide loop to keep the leg supports together. A small leg support opening minimises the risk of the user slipping through the sling.


A pocket in the lower back and on the leg supports facilitates application of the sling. The sling has a red seam in the middle of the back to help ensure that it is positioned symmetrically.

Anti-slip for added security

EvoSling Ampu is made of polyester net and has anti-slip material on the inside under the buttocks to provide added security during lifting. The user has a comfortable reclined sitting position. The leg supports are padded for optimum comfort.



Molift EvoSling Ampu MediumBack




Size

XS
S
M
L
XL
XXL

Net Padded Item no.:

1723210
1723220
1723230
1723240
1723250
1723260

| Size | Ⓐ Back height | Ⓑ Width |
|------|---------------|---------|
| XS | 47 cm | 21 cm |
| S | 52 cm | 26 cm |
| M | 63 cm | 32 cm |
| L | 77 cm | 42 cm |
| XL | 90 cm | 48 cm |
| XXL | 100 cm | 54 cm |





Molift EvoSling Comfort MediumBack

A FULL BODY SUPPORT SLING FOR USERS WHO NEED EXTRA SUPPORT

Support under the entire sitting area

EvoSling Comfort MediumBack Net is a full-body support sling without split leg support and with a shoulder-height back. The sling gives the user support under the entire sitting area and is ideal for users with joint problems or high sensitivity to pain.

The sling is positioned under the user in the lying position. It is made of polyester net, which provides good ventilation and can be left in the wheelchair underneath the user.

Adjustable sitting angle

Multiple hoisting loops on both the head and leg supports allow the sitting angle to be adjusted.


The sling has a red seam in the middle of the back to help ensure that it is positioned symmetrically.



Guide straps


The guide straps are built into the side straps to ensure that the sling is comfortable.

Molift EvoSling Comfort MediumBack



| Size | Net Item no.: |
|------|---------------|
| XS | 1722610 |
| S | 1722620 |
| M | 1722630 |
| L | 1722640 |
| XL | 1722650 |
| XXL | 1722660 |

| Size | (A) Back height | (B) Width |
|------|-----------------|-----------|
| XS | 77 cm | 85 cm |
| S | 90 cm | 86 cm |
| M | 102 cm | 94 cm |
| L | 115 cm | 113 cm |
| XL | 118 cm | 120 cm |
| XXL | 120 cm | 126 cm |



Molift EvoSling accessories

Molift EvoSling Support

Accessories for slings without head support



The head support is attached to the sling and its length can be adjusted as required using the plastic buckle. It gives the user comfortable support for the head and is a useful addition when the user requires more support.




Product
Molift EvoSling Support

Item number
1723400

Molift EvoSling Soft Padding

For pressure-sensitive users



EvoSling Soft Padding consists of pressure-relieving padded sleeves which are slipped over the leg parts of slings with split-leg support.

Product (sold in pairs)
Molift EvoSling Soft Padding
XS-M
Molift EvoSling Soft Padding
L-XXL
Molift EvoSling Hygiene Soft Padding
XS-XXL


Item number
1723710
1723810
1723740

Molift EvoSling Extension Loop

Straps for easier adjustment




Used to extend the hoisting straps or to keep the leg supports together on slings with split-leg support.



Product
Molift EvoSling Extension loop
16 (+9) cm
Molift EvoSling Extension loop
25 (+18) cm

Item number
1723510
1723520

Molift Strap Marker



For marking a selected hoisting loop for head or leg support hoisting straps. Mark the hoisting loops selected during individual testing of the sling.

Product (sold in pairs)
Etac Molift Strap Marker

Item number
1723600

Molift EvoSling FlexiStrap

Avoid heavy lifting



The Molift EvoSling FlexiStrap is used to easily and comfortably lift a leg or arm during care or handling. The hoisting strap is made of padded polyester and comes in two widths. There is a pocket on one side of the strap to facilitate positioning.



Product

Molift EvoSling FlexiStrap, standard
15 x 95 cm
Molift EvoSling FlexiStrap, wide
25 x 95 cm

Item number

1723310
1723330

Immedia SlingOn

When the user can't be turned onto their side




If you combine an "Immedia SlingOn" with a sling, you can position the sling directly under the user, i.e. without first having to turn the user onto their side.

Product

Immedia SlingOn, small
80 x 100 cm
Immedia SlingOn, medium
80 x 170 cm
Immedia SlingOn, large
120 x 150 cm

Item number

IM437/S
IM438
IM437





“Everybody is entitled to live a free and independent life and to fulfil their dreams
whatever their physical circumstances.”

You can also find news and more information at www.etac.com