

# EZ TRODE®

## Self-adhering, Reusable Electrodes

**Easy to Use**—Simply peel from protective backing and apply to the skin. After treatment is complete, return to the protective backing and pouch for storage.

**Excellent Conductivity**—Silver coating on electrode backing provides uniform current dispersion with high conductivity.


**Pliable**—Conforms to irregular body surfaces for complete contact, even during strenuous muscle stimulation protocols.

**Long Lasting**—Skin friendly conductive gel maintains smooth skin contact through many stimulation sessions.

**Foil-Lined Pouch**—Sealed pouch provides a moisture barrier to extend shelf life. Easy closing top keeps electrodes fresh after opening.

**New Soft Connectors**—Molded rubberized connectors assure easy insertion and removal of electrode cable-pins into the electrode connector.

**Available in Four Convenient Sizes**—EZ Trode electrodes come in four convenient sizes and shapes to provide maximum therapeutic flexibility.



- *2" Round*— (2221) Ideal for interferential stimulation of the ankle, knee, shoulder, neck and elbow.
- *2.75" Round*— (2222) A larger size electrode for interferential stimulation of the low back.
- *2" x 5" Rectangle*— (2223) Fits large muscles such as the quadriceps or hamstrings for muscle stimulation protocols.
- *2" Square*— (2224) Square electrode maximizes stimulation surface area while keeping the same price as 2" round electrodes. Rounded corners help reduce edge curl.