

# Multiglide Glove

## Short Instruction

**immedia**<sup>®</sup>

a part of Etac



**Material**

100% nylon, polyester

Class 1 - Medical device directive  
MDD 93/42/EEC

**Cleaning**

Do not add fabric conditioner  
Can be cleaned with a solvent-free cleaning agent with a pH level between 5-9, or a 70% disinfection solution.



**Item no**

**EAN no**






**Size in mm**

IM35  
IM8000701

5708012254576  
5708012266852

W200 x L500  
W150 x L350

# Multiglide Glove



## Check always:



## Contact / Download

Etac A/S  
[info@etac.dk](mailto:info@etac.dk)  
[www.etac.dk](http://www.etac.dk)

Etac GmbH  
[info@etac.de](mailto:info@etac.de)  
[www.etac.de](http://www.etac.de)

Etac A/S  
[hovedkontor.norge@etac.com](mailto:hovedkontor.norge@etac.com)  
[www.etac.com](http://www.etac.com)

R82 UK Ltd.  
[enquiries@etac.uk](mailto:enquiries@etac.uk)  
[www.etac.com/uk](http://www.etac.com/uk)

Etac AB  
[info@etac.se](mailto:info@etac.se)  
[www.etac.se](http://www.etac.se)

Etac Export  
[etac.immedia@etac.com](mailto:etac.immedia@etac.com)  
[www.etac.com/exp](http://www.etac.com/exp)

Etac R82 BV  
[info.holland@etac.com](mailto:info.holland@etac.com)  
[www.etac.com](http://www.etac.com)

R82 Inc.  
[information@r82.com](mailto:information@r82.com)  
[www.r82.com](http://www.r82.com)

## Manufacture

Etac Immedia A/S  
 Egeskovvej 12  
 DK-8700 Horsens  
[etac.immedia@etac.com](mailto:etac.immedia@etac.com)

 **etac**<sup>®</sup>  
 Creating possibilities