

Patient Lift Sling

Residents who need assistance with mobility will benefit from the lift sling, which is compatible with other manufacturer's lifts. With a variety of styles and models, caregivers can select the type that best suits their resident's needs.

Mesh Full Body Sling

Full body slings are ideal for residents who need full support, especially for the head and neck. The mesh material offers more breathability and dries easily when exposed to moisture.

Mesh Full Body Sling with Commode

Full body sling with commode option provides easy toileting access, and can also be used for residents who need to be in a constant seated position.



Solid Full Body Sling

Polyester full body slings are more durable and provide full body support. This sling is ideal for transfer from a wheelchair to a bed or another chair, or transfer from floor to bed.

Patient Lift Sling

Solid Divided Leg Sling

The divided leg sling is the most versatile option, and is quick and easy to use. It helps the residents distribute their weight more comfortably.



Sit to Stand Sling

The Sit to Stand sling is used for transfer from a seated position in and out of bed to a chair, wheelchair, toilet or from chair to chair. It is designed for residents with minimal weight-bearing ability, who have head and neck control, can bend at the hips, knees and ankles, and can sit up on the edge of the bed.

Model #	Type	Size	Binding Color	Dimension A (in)	Dimension B (in)	Dimension C (in)	Weight Capacity
#71256	Full Body Sling with Commode, Mesh	M	Purple	39	51	24	600 lbs
#71257		L	Green	43	55	24	600 lbs
#71258		XL	Blue	43	59	27	600 lbs
#71259		XXL	Black	43	51	44	600 lbs

