

ProMotion® Plus

Orthotic Insoles

Full Length PORON® Cushioning Insoles provide relief from:

- › Over-pronation
- › Plantar fasciitis/heel pain
- › Arch pain
- › Ankle and knee pain



DARCO

DARCO ProMotion® Plus Orthotic Insoles:

- › The soft high filament count polyester top fabric contains a permanent anti-microbial for odor control and low friction, allowing freedom of movement without socks bunching up.
- › Full length PORON® Cushioning foam layer provides superior shock absorption.
- › Heel to toe cushioning foam that won't flatten or bottom out. Superior feel and long term wear benefits.
- › EVA covered polypropylene heel and arch support. Provides a firm support that cradles the heel and arch area, The EVA surface prevents the insole from shifting while walking and adds durability.

Indications: Over-pronation, plantar fasciitis/heel pain, foot and arch pain, ankle and knee pain.

Preferred Shoe Type: Dress, casual, and athletic shoes.


Directions: Remove existing insoles from shoes and replace with Darco ProMotion Plus Insoles.

Care Instructions: Wipe with a damp cloth. Do not submerge in water. Allow insoles to air dry.

Warning: Always consult a medical professional if persistent pain, discomfort or swelling continues.

Not intended for severe foot deformities or circulation problems, infections or foot ulcers.

Suggested Code: L3040



PORON is a licensed trademark of Rogers Corporation.