



CORRXIT®

with Ambulatory attachment



Washable	Measure Heel to Great Toe	Maximum Calf Circ.	Model #
Average	8 1/2 – 9 3/4" (21.6 – 24.8 cm)	17" (43.2 cm)	03CX-A
Large	9 3/4 – 11" (24.8 – 27.9 cm)	17" (43.2 cm)	03CX-L
Bariatric	9 3/4 – 11" (24.8 – 27.9 cm)	20" (50.8 cm)	03CX-B

Assigned PDAC Code L4396 or OTS L4397

- The unique heel design of the splint eliminates pressure points and allows for visual examination and circulation of air over the heel.
- The Toe Guard reduces bed linen skin irritation by separating the bedding from contact with the foot.
- The Stabilizer Bar effectively prevents rolling of the patient's leg and controls internal or external rotation
- The washable, anti-microbial foam liner of the CORRXIT® wicks away moisture, keeping the skin dry and comfortable
- The low profile Ambulatory Attachment (shown) is easily applied and removed without the use of tools