

THERMAL CARPAL TUNNEL GLOVE

- Commonly used in the treatment of carpal tunnel syndrome
- Reduces wrist motions that may cause carpal tunnel syndrome
- Supports the wrist while easing arthritis pain
- Textured palm and fingers for non-slip
- Open fingertips for improved sense of touch
- Contoured aluminum support significantly reduces wrist motion
- Adjustable closure for comfort

Swede-O Thermal Carpal Tunnel Glove is designed to provide support for weak, injured, or aching wrists while helping to reduce symptoms of arthritis,



poor circulation and other hand conditions. The thermal lining captures and retains your natural body heat to help increase circulation and aid in long lasting pain relief.

Available From:

FUNCTIONAL & COMFORTABLE

The Swede-O® Thermal Carpal Tunnel Glove provides comfortable support to the wrist and is commonly used in the treatment of carpal tunnel syndrome. The removable palmar spoon is preformed with a 25° angle to support the wrist in a neutral position.

The thermal lining captures and retains your natural body heat to help increase circulation and aid in long lasting pain relief. This form fitting glove provides mild compression to help decrease swelling and pain of the joints. The non-slip, textured palm and fingers with open fingertips are convenient for everyday activities. Specify left or right hand when ordering.



Materials: Outer lining: Nylon, Mid Layer: Polyurethane, Inner Lining:

Nylon & Spandex.

Latex: This product is not made with

natural rubber latex.

Product Selection:

SPECIFICATIONS

WST-6839 Swede-O® Thermal Carpal Tunnel Glove

Product Weight: .25 lb Dimensions: 6" x 9" x 2"

Case Weight: 3 lbs

Case Dimensions: 12" x 10" x 8"

LEFT	RIGHT	SIZE	MEASUREMENT	
WST-6839-L-1XS	WST-6839-R-1XS	XS	6-6 3/4"	15-17cm
WST-6839-L-SML	WST-6839-R-SML	S	7-7 3/4"	18-20cm
WST-6839-L-MED	WST-6839-R-MED	М	8-9"	20-23cm
WST-6839-L-LRG	WST-6839-R-LRG	L	9 1/4-10 1/2"	24-26cm
WST-6839-L-1XL	WST-6839-R-1XL	XL	10 3/4-11 1/2"	27-29cm
WST-6839-L-2XL	WST-6839-R-2XL	2XL	11 3/4+	30cm +

Your Thermal Support should be a snug and comfortable fit. Not too tight, not too loose. If in doubt, try next larger size.

Measure the circumference of the hand around the knuckles

