


MHRTSL

Raised Toilet Seat with Legs Assembly Instructions





MHRTSL

Raised Toilet Seat with Legs

Assembly Instructions

1. Push the four legs into the seat
2. Insert the two safety bars into the four holes in the top of the plastic seat.
3. Adjust the legs so the seat sits firmly on the bowl and the legs rest on the floor.
4. The seat is now ready for use. It may be cleaned by using a mild detergent solution and sterilized with alcohol.

Note: Screws can be used as an option. Use in four holes.



MOBB Health Care Ltd.
116A Viceroy Rd. Concord,
ON L4K 2M1 Canada