

Walkers

Decreased mobility doesn't have to mean the loss of independence. McKesson walkers are a convenient and affordable mobility solution for users who have difficulty with balance or standing, suffer arthritic pain, or experience fatigue or shortness of breath from walking. Walkers are commonly preferred for stability and balance. McKesson walkers feature lightweight construction and a push-button release mechanism for easy folding and storage. Choose McKesson walkers for durability, convenience, and comfort.

Features and Benefits

- Sturdy, 1-inch, aluminum A-frame construction is lightweight without sacrificing durability for patients up to 300 lbs
- Steel legs and side braces maximize durability and strength to support patients up to 500 lbs (applicable models only)
- Push-button, dual-release mechanisms for easy-folding using fingers, palms or side of hand to accommodate users' abilities and travel convenience
- Rubber-like tips made of polyurethane blend to provide better grip than standard walkers
- 5-inch front wheels and rear glide caps allow walker to slide easily and smoothly over most surfaces (only for walkers with wheels)
- Adjustable height and vinyl-contoured hand grips
- U-shaped frame provides greater clearance and the cross frame design offers improved stability
- Available in adult and junior sizes
- Not made with natural rubber latex



The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.