

## Airex® Mats

Airex® Mats are the finest lightweight, portable mats available. Firm, comfortable, and colorful, they have excellent cushioning qualities. Being closed cell, these mats are impervious to water and are easy to clean with a little soap and water. Traverse ribbed surface pattern resists tearing and allows the mat to be rolled up easily. Airex® Mats readily lie flat and firmly grip the floor or tabletop to provide a stable, nonslip surface. Tiny air-filled cells act as individual air cushions that firmly resist the “bottoming out” of a knee or elbow pushed into the mat. These mats are impervious to water and are easy to clean. All mats are sanitized; a molded in process that guarantees continued hygiene. They do not contain latex.



**Mat Bags** are made of lightweight nylon with wrap-around zippers and tough web handles. Simple and easy to use, **Mat Straps** with Velcro® fasteners hold a rolled mat neatly for storage or carrying. Use one strap for narrow mats, two for wide mats.

TOP TO BOTTOM: ATLAS, HERCULES, CORONA AND CORONELLA MATS

### Airex® Mats, Choose Color

THA-MAT-A-G Atlas Mat – Green  
THA-MAT-A-R Atlas Mat – Red

Dimensions: 48"W x 78"L x 5/8"H

THA-MAT-B-C Mat Bag for Corona  
THA -MAT-B-CL Mat Bag for Coronella  
THA-C-B Corona Mat – Blue  
THA-C-G Corona Mat – Green  
THA-C-R Corona Mat – Red

Dimensions: 39"W x 72"L x 5/8"H

THA-MAT-CL-B Coronella Mat – Blue  
THA-MAT-CL-G Coronella Mat – Green  
THA-MAT-CL-R Coronella Mat – Red

Dimensions: 23"W x 72"L x 5/8"H

THA-MAT-H-B Hercules Mat – Blue  
THA-MAT-H-G Hercules Mat – Green

Dimensions: 39"W x 78"L x 1"H

THA-MAT-S-M Mat Strap – Medium

### Hop 'N Pop®

THA-HNP-100 Hop 'N Pop® Mat

Overall dimensions: 2 1/2 ft. W x 8 1/2 ft. L

THA-HNP-100-R Replacement Bubble Sheets (set of 6)



Hop 'N Pop® is a sensory-based, large exercise mat that encourages active play for pediatric adaptive physical and occupational therapy programs. The activity fitness mat contains a brightly-colored, custom-sized sheet of close to 2,000 bubbles that snap and pop each time a child jumps, hops, dances, or plys on the mat, thus helping children discover the relationship between cause and effect. The mat's top layer is made of a sturdy, washable, mesh fabric with vibrant graphic designs, while the bottom layer is a 1/4" thick, non-skid foam vinyl, which is heavy duty for extra stability and absorbing impact.

Each complete set comes with four rolls of custom-sized bubble sheets, an adjustable shoulder strap for portability, a Hop 'N Pop® key chain, two keys to the child-proof zipper lock, instruction and care sheet, game booklet, and warranty. Replacement **Bubble Sheets** are available in sets of six.



### Sensory Tactile Mat™

THA-KIN-800TM Sensory Tactile Mat™

Dimensions: 18 1/2"W x 18 1/2"L x 1 1/2"H (each piece)

Age range: 2 yrs. and up

## Sensory Tactile Mat™

A lively and vividly shaped mat may be pieced together and the different grain designs enhance tactile sensitivity. Multicolored, with 6 interlocking squares in a set.