

THE ASCEND & GO SIT TO STAND LIFT

This transfer lift is ideal for transferring patients with difficulty walking, or patients who are a fall risk. It encourages early mobility movement by encouraging patients to be active during sitting to standing transfers.

Best Option for Partial Mobility



400 LBS.
WEIGHT CAPACITY



50 LBS.
UNIT WEIGHT



ADJUSTABLE
KNEE SUPPORT



◆ Pivoting Seat

Makes settling into the lift easy and secure, allowing patients to remain in their preferred standing position during transfer.

◆ Knee Supports

Provide extra stability, security, and comfortability while being transferred in the seated position.

◆ Dual Locking Wheels

Designed to ensure safety and unwanted movement during transfer into the device, or during transit.

PRODUCT SPECIFICATIONS

SKU: MP-P02

LIFT CAPACITY	400 LBS.
UNIT WEIGHT	50 LBS.
PLATFORM DIMENSIONS	19" LENGTH x 19" WIDTH
ASSEMBLED UNIT DIMENSIONS	32.2" LENGTH x 20.8" WIDTH x 41.7" HEIGHT
KNEE PAD MEASUREMENTS	15.3" WIDTH x 10.2" HEIGHT
CHASSIS EXTENSION	19.2" - 29.9"