

PO: _____	Date: _____	Customer #: _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



Standard Features	
<ul style="list-style-type: none"> · Tubular Steel, Powder Coated Frame · 45 Degrees, Infinitely Adjustable Seat Tilt · Fixed Height Armrests · Adjustable Seat Depth 16" - 20" · Adjustable Seat to Floor Height 14" - 19" · Choice of Wheel, Caster and Fork Sizes · Push Lock Brakes 	<ul style="list-style-type: none"> · Drop In, Swing Away Foot Rests · Composite Foot Plates · Padded Side Panels · Adjustable Backrest Angle · 250 lbs. Weight Capacity · Seat cushion, back and headrest not included

Comfort Rehab Wheelchair (587) - E1161

HCPCS codes supplied are for reference only and do not ensure coverage or payment. The HCPCS provided are based on our interpretation of the code definitions. HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

Chair Color

Frame Color: Black Blue Bronze Nickel Orange Red Sand Yellow

Seat and Back Options

Seat Width: 16" 18" 20" (E2201) 22" (E2201)

Seat Depth: 16" 17" 18" 19" 20" (E2203)

Solid Seat Pan (E2231) Factory Installed

Omit Seat Pan

Seat to Floor Height: 14" 14.5" 15" 15.5" 16" 16.5" 17" 17.5" 18" 19"

Back Cane Height: Ultra-Low Double Bend (16") Fixed

Ultra Low (16") Fold Down


Low (18") Fold Down

Medium (21") Fold Down

Tall (24") Fold Down

Back Cane Angle: 5° 10° 15° 20°

Refer to Seat Height & Wheel Combinations Matrix for assistance.



Front Mount Options	Wheel Options	Rear Mount Options
<input type="radio"/> 5" Fixed Front Casters (pair)	<input type="radio"/> 20" Rear Mounted Mag Wheels (pair)	
<input type="radio"/> 6" Fixed Front Casters (pair)	<input type="radio"/> 22" Rear Mounted Mag Wheels (pair)	
<input type="radio"/> 8" Fixed Front Casters (pair)	<input type="radio"/> 24" Rear Mounted Mag Wheels (pair)	

Caster Forks: Required

Small - Recommended for Small Mag Wheel (20" or 22") & Lower Seat to Floor Heights

Medium - Recommended for Large Mag Wheel (24") & Higher Seat to Floor Heights

Refer to Seat Height & Wheel Combinations Matrix for assistance.

Mag Wheel Lock Brake Extension Handles (pair) (E0961)

Anti-Tippers (E0971) - Required Safety Features - Factory Installed

Armrest Options

Fixed Height Armrests

Height Adjustable, Removable Armrests (E0973)

Hanger & Legrest Options	Legrest Accessories:
Legrest Hanger Style: must choose one	<input type="radio"/> Heel Loops (pair) E0951
<input type="radio"/> Drop In Straight (pair)	<input type="radio"/> Legrest Panel (for calf support) E0995
<input type="radio"/> Drop In Flared (pair)	<input type="radio"/> Foot Board
<input type="radio"/> Pin and Saddle Straight (pair)	<input type="radio"/> F-SW-ELV Swing Away, Elevating Legrests (E0990) (pair)
<input type="radio"/> Drop In IKON Euro Style (pair)	<input type="radio"/> F-SW-ELV-IKN Swing Away, Elevating - Euro Style Legrests (E0990) (pair)
<input type="radio"/> F-SW75ADJ-587V3 Swing Away 75° Footrests (pair)	<input type="radio"/> w/ Angle Adjustable Footplates (K0040) (pair)
<input type="radio"/> F-SW60ADJ-587V3 Swing Away 60° Footrests (pair)	
<input type="radio"/> F-SW-IKN Swing Away Footrests (silver) (pair)	
<input type="radio"/> F-SW-IKN-ADJ Swing Away Footrests (Black)	
<input type="radio"/> w/ Angle Adjustable Footplates (pair) (K0040)	

Common Accessories

Pelvic Strap: Pelvic Positioning Strap Grommet Mount (E0978) Ctr Buckle Side Buckle

Trays: ABS Tray with Strap (E0950) (16", 18" 20" 22" - Matches Corresponding Seat Width) Sewn Strap Grey Black

Padded ABS Tray Cover (ABS Tray only)

Clear Tray with Strap (E0950) (16", 18", 20" - Matches Corresponding Seat Width) Sewn Strap

Other: IV Pole & Holder (includes mounting bracket) - (K0105) Front Mounted OR Rear Mounted

O2 Holder (includes mounting bracket) (E2208)

Accessory Bag (Black) Each

Accessory Bag Personalizing 2 Initials

Cup Holder

Brush Font OR Block Font

Matrx® Seating Series



Matrx® Seat Cushions

For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions

Seat Cushion Depth: 16" 18" 20" 22" (Only available w/ PS Cushion)

Seat Cushion Width: 16" 18" 20"

Libra (E2624) - Recommended for highest risk of or existing skin breakdown
 Vi (E2607) - Recommended for moderate to high risk of skin breakdown
 PSP (E2607) - Recommended for moderate to high risk of skin breakdown
 PS (E2605) - Recommended for low to moderate risk of skin breakdown

Seat Cushion Width: 22"

Libra (E2625) - Recommended for highest risk of or existing skin breakdown
 Vi (E2608) - Recommended for moderate to high risk of skin breakdown
 PS (E2606) - Recommended for low to moderate risk of skin breakdown

Libra Cushion Accessories: **Libra Fluid Supplement**
 Libra 1" Hemi Wedge
 Libra 1" Anterior Seat Wedge
 Libra Tapered Wedge Left
 Libra Tapered Wedge Right
 Libra Abductor

Vi, PSP & PS Cushion Accessories: **Pelvic Obliquity Kit**

Cushion Rigidizer:	E0992	16" Depth	18" Depth	20" / 22" Depth
All Cushions		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
		N/A	<input type="radio"/>	<input type="radio"/>

Matrx® Backs

For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions

Back Height: 16" 18" (excludes PB & PB Deep) 20"

Back Width: 16" 18" 20"

Elite TR (E2620) - Recommended for tilt in space & recline applications
 Elite (E2620) - Recommended for more active users
 Elite Deep (E2620) - Recommended for active users with contoured support needs
 PB (E2615) - Recommended for less active users seeking comfort
 PB Deep (E2615) - Recommended for less active users seeking comfort & support

Back Width: 22"

Elite TR (E2621) - Recommended for tilt in space & recline applications
 Elite (E2621) - Recommended for more active users
 Elite Deep (E2621) - Recommended for active users with contoured support needs

Lateral Support Accessories: **Elite Back Fixed Laterals 5"Hx7.5"L (E0956)**
 Elite Back Swing away Lateral Hardware (E0128)
 Swing away Lateral Pad (each) 4.25"Hx6.5"L (E0956)

Privacy Flap Accessories: **16" Width Small Privacy Flap**
 18" Width Medium Privacy Flap
 20" Width Large Privacy Flap
 22" Width Extra Large Privacy Flap

Matrx® Headrest

Headrest Mounting: **Elan Headrest Hardware (E0128)** - Required with Selection of Headrest

Headrest Options: **Elan Standard Pad with Standard Cover (E0955)** 10" 14"
 Elan 4 Point Pad (E0955) 9" x 7" 11" x 10"
 Elan Headband with 4 Point Pad (E0955) Small Large

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$ _____

587V3: SEAT HEIGHT & WHEEL COMBINATIONS									
SEAT HEIGHT	MAG WHEEL	CASTER WHEEL	CASTER FORK	MAG PLATE HOLE #	CASTER FORK HOLE #	ANTI-TIPPER MOUNT	ANTI-TIPPER HOLE #	SEAT TAB	
14"	20	5	SMALL	1	1	TOP	3	2	
	22	5	SMALL	1	2	TOP	2	1	
15.5"	20	6	MEDIUM	2	1	TOP	2	1	
	22	6	MEDIUM	1	1	TOP	2	1	
	24	6	SMALL	1	2	TOP	1	1	
	24	5	MEDIUM	1	2	TOP	1	1	
16"	20	5	SMALL	3	3	TOP	1	1	
	20	5	MEDIUM	3	3	TOP	1	1	
	20	6	MEDIUM	3	2	TOP	1	1	
	22	5	SMALL	2	3	TOP	1	1	
	22	6	SMALL	2	2	TOP	1	1	
	22	5	MEDIUM	2	3	TOP	1	1	
16.5"	20	6	SMALL	3	3	TOP	1	1	
	24	6	SMALL	2	3	BOTTOM	3	1	
	24	5	MEDIUM	2	4	BOTTOM	3	1	
	24	6	MEDIUM	2	3	BOTTOM	3	1	
17"	20	5	MEDIUM	4	5	BOTTOM	3	1	
	20	6	MEDIUM	4	4	BOTTOM	3	1	
	22	5	MEDIUM	3	5	BOTTOM	3	1	
17.5"	22	6	MEDIUM	3	4	BOTTOM	3	1	
	20	8	SMALL	4	3	BOTTOM	3	1	
	22	8	SMALL	3	3	BOTTOM	3	1	
	24	8	SMALL	3	3	BOTTOM	2	1	
18"	24	6	MEDIUM	3	5	BOTTOM	2	1	
	24	8	MEDIUM	3	3	BOTTOM	2	1	
	20	8	MEDIUM	5	4	BOTTOM	2	1	
19"	22	8	MEDIUM	4	4	BOTTOM	2	1	
	24	8	MEDIUM	4	5	BOTTOM	1	1	

Hole numbering begins at the top of the plate, fork and anti-tipper. Hole number 1 is the uppermost hole.
 Do not use the rear most vertical row/column of holes on the mag plate for 20" mag wheels

Form08CONFIG-404 587V3 Rev. 1