

Deluxe Clinical Care Recliner

The Deluxe Clinical Care Recliner sets a new standard of comfort for the patient and ease of use for the caregiver.

- Available With Heat and Massage Option.
- The Series FR566G and FR566GH Deluxe Clinical Care Recliners feature three positions for patient comfort. The fourth position is a Trendelenburg (or Shock Position) that must be activated by a caregiver. The Trendelenburg Position cannot be activated by a patient while in the recliner.
 - The first position is the **Upright Position**. This position is ideal for patient transport and can be used to initiate patient treatment / therapy.
 - The second position is the **First Recline Position** and is also known as the “TV” position. This position is ideal for patient treatment, rest and reading.
 - The third position is the **Second Recline Position** and is also known as the “Full Recline” position. This position is ideal for patient treatment, rest and reading.
- The wide, contoured seat and back, along with an integrated footrest, provide an enhanced level of comfort. Built-in arm pads in the back give improved support for the arms.
- Wall saver armrest caps increase durability of the armrests and protect walls from marring.
- Higher floor to back height enables the caregiver to raise the patient up with less effort, reducing the chance for a back injury.
- Extra deep Trendelenburg Position is easily actuated from either side of recliner by a conveniently located foot pedal. Steeper angle helps improve patient stabilization.
- Four Tente® swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with easy to locate black locking tab lock wheel and swivel with single motion. Directional wheel lock with high visibility green tab allows recliner to track straighter with less effort during transport.
- Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Easy to actuate by caregiver. Side table surface is smooth and features tapered edges for easy cleaning.
- Extra heavy duty recliner mechanism.
- 350 lb maximum weight capacity, EVENLY DISTRIBUTED.



Upholstered inside panels

Smooth tables with no-seam tapered edge

Specifications

Overall Height.....	50"	Seat Depth.....	20"
Overall Width.....	31"	Depth Fully Reclined	
Overall Depth.....	43"	In Trendelenburg	75"
Width Between Arms	21"	Width With Both Side Tables	
Arm Height From Floor	27.5"	in Upright Position	51"
Arm Height From Seat.....	8"	Shipping Weight	145 lb
Seat Height.....	22"	Max. Weight Capacity.....	350 lb
Seat Width.....	20.5"		



Easily removed side panels enable easy access to recliner for cleaning and maintenance.



Easy-to-achieve Trendelenburg position. Trendelenburg position can be actuated from either side of the recliner by a convenient foot pedal.



Multiple Positions for Enhanced Patient Comfort



Colors may vary slightly from actual vinyl

Standard	with Heat & Massage	
FR566G401	FR566GH401	Vintage Gold
FR566G409	FR566GH409	Warm Taupe
FR566G427	FR566GH427	Blue Ridge
FR566G432	FR566GH432	Imperial Blue
FR566G454	FR566GH454	Royal Blue
FR566G6705	FR566GH6705	Dolce Sand
FR566G6714	FR566GH6714	Dolce Burgundy
FR566G6717	FR566GH6717	Dolce Moss
FR566G6724	FR566GH6724	Dolce Jet
FR566G6725	FR566GH6725	Dolce Indigo
FR566G6727	FR566GH6727	Dolce Meteor
FR566G8504	FR566GH8504	Sand
FR566G8516	FR566GH8516	Chestnut
FR566G8533	FR566GH8533	Port
FR566G8561	FR566GH8561	Mulberry
FR566G8567	FR566GH8567	Dove
FR566G8584	FR566GH8584	Taupe
FR566G8602	FR566GH8602	Oatmeal
FR566G8604	FR566GH8604	Concrete
FR566G8611	FR566GH8611	Terra Cotta
FR566G8615	FR566GH8615	Berry

Standard	with Heat & Massage	
FR566G8622	FR566GH8622	Grey Moss
FR566G8624	FR566GH8624	Granny Smith
FR566G8626	FR566GH8626	Eucalyptus
FR566G8628	FR566GH8628	Sea
FR566G8801	FR566GH8801	Marine Blue
FR566G8805	FR566GH8805	Jet
FR566G8808	FR566GH8808	Champagne
FR566G8809	FR566GH8809	Taupe
FR566G8818	FR566GH8818	Sage
FR566G9201	FR566GH9201	Cobblestone
FR566G9202	FR566GH9202	Periwinkle
FR566G9205	FR566GH9205	Gypsum
FR566G9207	FR566GH9207	Ice Blue
FR566G9208	FR566GH9208	Wineberry
FR566G9210	FR566GH9210	Cypress
FR566G9212	FR566GH9212	Blue Jade
FR566G9214	FR566GH9214	Steel Blue
FR566G9715	FR566GH9715	Eucalyptus
FR566G9722	FR566GH9722	Desert Sage
FR566G9764	FR566GH9764	Steel Grey