

DORSAL

Night Splint



 DeRoyal®

Finally, a splint the patient will wear.

Help your patients get back on their feet with this

Dorsal Night Splint

- Thermoplastic design holds foot in 90° neutral position during sleep to help alleviate morning pain and discomfort associated with plantar fasciitis
- Plastic free bottom and non-skid transfer surface allows for safe minimal ambulation when necessary
- Neoprene straps with hook and loop closure around foot and calf
- Low profile for cooler, more comfortable fit
- Universal right/left

Indications

Plantar fasciitis, Achilles tendonitis,
Tendonosis, Drop foot, Post-static pain



Ordering Information

| Product # | Size | Male Shoe Size | Female Shoe Size |
|-----------|---------------|----------------|------------------|
| NS2000-20 | Small/Medium | 4 - 9 | 5 - 10 |
| NS2000-21 | Large/X-Large | 9½ - 14 | 10½ - 15 |

**Product is not designed for walking outdoors or on slippery, uneven surfaces*

Diverse. Dynamic.

