

The Duet has locking flip up armrests and a flip down footrest to aid the user while the Duet is used in the transport chair mode. A safety lock release is used to store the arms and a velcro strap is used to store the footrest while using the Duet for ambulation.

The tool-free two position contoured padded back rest is reversible depending on whether you are using the Duet in an ambulatory or transport chair mode.

- Easy release loop locks.
- Comfortable padded seat with a large underseat pouch for safe storage of personal items. Inside the pouch is a handle to aid with the folding.
- The folding frame has a lock out feature to insure the unit stays open while in use and insures the unit will stay folded during transport.
- Large 8 inch casters are used for both indoor and outdoor activities.
- Weighs only 19 lbs.

SPECIFICATIONS

Limited Lifetime Warranty.

SEE PAGE 52 & 53 FOR REPLACEMENT SEATS.









EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB.



Blue Burgundy Black

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Item #	Overall Width (open)	Overall Height	Seat To-Floor Height	Handle Height Range	Back Height	Width Seat	Depth Seat	Overall Length	Product Weight	Carton Shipping Weight	Weight Limit
DUET	25.5"	37"	21 "	31.5"- 37"	36"	13.5"	13"	26"	19 lbs.	23 lbs.	300 lbs.