

ThermalSoft Durapak hot and cold packs maintain an even temperature over the entire pad. Their soft gel contents are leakproof and contain no liquids or dyes. ThermalSoft Durapak hot and cold packs are biodegradable and non-caustic. And, because ThermalSoft Durapaks are economical and reusable, you can send them home with your patients for continued therapy at home between office visits.

For heat therapy, microwave or heat the pack in water. ThermalSoft Durapaks keep their pad like shape even at high temperatures and will not liquefy when heated.

When used for cold therapy, ThermalSoft Durapaks stay flexible when frozen. They bend to fit around treatment areas and maintain their thickness no matter how severe the bend.

## **Durapak Features**

- Leakproof
- Economical
- Use for **Both** hot and cold
- Freeze for cold
- Microwave for heat
- No need to heat in a bowl of water
- Soft and flexible





## **Durapak Sizes**

- 9841—Small Pack (4" x 6")—48 per case
- 9842—Half Size Pack (5" x 10")—24 per case
- 9843—Back Pack (8" x 11")—24 per case
- 9844—Extra Large Pack (12" x 15")—12 per case