

Encore Rehab Wheelchair Model: 500RW

	•				~
PO:		Date:			Customer #:
Customer Nam	e:			Ship To Name:	
Address:				Address:	-
Phone:	-			Phone:	
Contact:				Contact:	
Email:		_		Tag:	-
				Standard	Features
		Solid Seat Pan			• 250lb Weight Capacity
		Tubular Steel, Powder Coated Frame			Adjustable Seat to Floor Height
		• Seat cushion, back & headrest not included			• 20" Front Mount Mag Wheels
	Discould be a second	 Independent Height Adjustable Arms Rocking Mechanism (allows up to 33° of tilt) 			 2 x 5" Rear Total Lock Casters Seat Mount Hangers
	C C C C C C C C C C C C C C C C C C C	• Fold Down and Angle Adjustable Back Canes			Removable Arms
	and in	• 33° Infinitely Adjustable Seat Tilt			Tomo radio : Imp
		Encore Rehab Wheelchair® -			1161
C 4 XXV 141	O 168 O 108	Seat ar	nd Back	•	0.168
Seat Width:	○ 16" ○ 18"			Seat Depth:	○ 16" ○ 17"
	20" (E2201) 22" (E2201)				○ 18" ○ 19" ○ 20" (E2203)
	○ 22 (E2201)			O Solid Seat Pa	<u> </u>
			•	Omit Seat Pa	
				Back Cane	Ultra-Low Double Bend (16") Fixed
Seat Height:	○ 14"			Height:	O Low (18") Fold Down
~ ····	<u>)</u> 15"			g	Medium (21") Fold Down
	○16" ○17" ○18"	○19"			◯ Tall (24") Fold Down
	Front Options	Whe	eel Optio	ns	Rear Options
○ 20" Front	Mounted Mag Wheels			∫ 5" Total Lock	k Rear Casters
	Mounted Mag Wheels				n Non-Locking Rear Casters
Front Ant	i-Tippers (E0971) Requir	•			Extension Handles (pair) (E0961)
		Arm	nrest Opt		
Armrest	O Soft Armrest			Armtop	Height Adjustable Arms (E0973)
Options:	Arm Trough (each)	:-1		Style:	*Height and Depth Adjustable (pr)
	O Padded Armrests (pa	ir)		Armtop Mount	*Not compatible with 22" mag wheels in the 14", 15", & 16" seat heights :: Seat Mounted Arms
		C	hair Col		Scat Wounted Airns
Frame Color	: O Black Orange	Carmine Red	mun con	Tungsten Sar	nd Ocean Blue
	○ Bronze ○ Satin N			<u> </u>	
	Drop In Style		Choose O	ne Hanger Styl	e Pin and Saddle Style
○ SM-HANG	GER-48RV4-RKR-IKN (p	air)		○ SM-HANG	GER-48RV4-RKR-PIN (pair) Less \$54
	Drop In Style	Legr	est Optio		Pin and Saddle Style
	Swing Away Silver (pair	r)		○ F-SW Swi	ing Away Footrests Black (pair)
	ADJ Swing Away Black				V Elevating Swing
	djustable Footplates (pair)				otrests Black(pair) (E0990)
	IKN Elevating Swing Away				OT-ENCORE Leg/Foot Board (E0995)
	djustable Footplates (pair) EL-ENCORE Legrest Pand			U LEGPAN.	EL-ENCORE Legrest Panel (E0995)
_	e	` '			
*Slight ground	interference may occur in the	lowest seat height. Slightly tilt chair or elevate legre	est as neede	ed	
		Comm	non Acce	ssories	
Seat Belts:	○ S-BELT	Pelvic Positioning Strap (E0978)			○ Ctr Buckle ○ Side Buckle
	○ S-BELT-A	Pelvic Positioning Strap Auto Buckle (E097)			○ Ctr Buckle ○ Side Buckle
T	○ S-BELT-587V3	Pelvic Positioning Strap Grommet Mount (E	E0978)		Octr Buckle Oside Buckle
Trays:	○ TRAY-ABS	Tray ABS with Strap (E0950)	050)		Sewn Strap Grey Black
	○ TRAY-ABS-OVR	Tray ABS with Strap (24" - 28" width) (E09	93U)		Sour Troy Strop
	○ TRAY-CLEAR	Tray Clear with Strap (E0950) Half Tray (E0950) (side specific)(with upper	r sida na	<i>1</i>)	Sewn Tray Strap
	○TRAY-HALF ○TRAY-COVER	Tray Cover (padded)(for TRAY-ABS)	г мае рас	<i>i j</i>	◯ Grey ◯ Black
Other:	() IV-POLE (K0105)	IV Pole & Holder (includes mounting bracket	et) - (F	ront or ORea	ar Mount
omer.	O2-HOLDER (E2208)	O2 Holder (includes mounting bracket) -			
	○ STABIL-FR	Rear Stabilizer Frame	, 11011t C	<u> </u>	
			Padding	Upgrades	
Positioning:	○ SR	Side Reduction Pads (pair) (E0956)			

Encore Rehab Wheelchair Model: 500RW

			Matrx	® Seating Series			-matr				
			Mada	® Seat Cushions			—Шил	' B			
		For Correct		B Seat Cusnions hould Match Selected Chair S	ant Dimonnia	***	SEATING SERIES				
Good Good Cooper	<u>()</u> 16"		20"	_		ns					
Seat Cushion Depth:	<u> </u>	<u>()</u> 18"	<u> </u>	22" (Only available w/ PS Co	usilioii)						
Seat Cushion Width:	<u>()</u> 16"	<u> </u>	○ 20"	Seat Cushion Widtl	h:	O 22"					
Libra (E2624)-Recommended				Libra (E2625) - Recommended for highest risk of or existing skin breakdow							
Vi (E2607) - Recommended for				Vi (E2608) - Recommended for moderate to high risk of skin breakdown							
O PSP (E2607) - Recommended			/I	O. P.G. (774-64-6)							
PS (E2605) - Recommended f	for low to moderate	risk of skin breakdown		OPS (E2606) - Reco	mmended for lov	w to moderate	risk of skin breakdown				
				Vi, PSP & PS	Cushion						
Libra Cushion Libra Flu	Accessor		O Pelvic O	Obliquity Kit							
Cusinon	Hemi Wedge										
◯ Libra 1"	Anterior Seat W	/edge		Cushion Rigidizer:	E0992	16" Depth	18" Depth	20" / 22" Depth			
O Libra Ta	pered Wedge Le	eft		All Cushions	16" Width	0	0	0			
◯ Libra Taj	pered Wedge Ri	ght			18" Width	0	0	0			
O Libra Ab	oductor				20" Width	0	0	0			
					22" Width	N/A	0	0			
			3.5	4 A.D. I							
		Ean Comp		atrx® Backs ould Match Selected Chair Sec	at Dim anai an						
Back Height: 016"	18"(excludes		20"	ula Maich Selectea Chair Sel	n Dimensions	•					
Back Height. 010	O 10 (excinues)	ть ать всер)									
Back Width: 016"	<u> </u>	○ 20"		Back Width: () 22"							
Elite TR (E2620) - Recommen	nded for tilt in spac	ce & recline applications		○ Elite TR (E2621) -	Recommended	for tilt in spac	e & recline applications				
Elite (E2620) - Recommended	d for more active us	sers			 ○ Elite TR (E2621) - Recommended for tilt in space & recline applications ○ Elite (E2621) - Recommended for more active users 						
Elite Deep (E2620) - Recomm			ort needs	 Elite Deep (E2621) 	- Recommende	d for active u	sers with contoured suppo				
		s seeking comfort						Ol'			
PB (E2615) - Recommended to PB Deep (F2615) - Recommended	nded for less active		> support					DF			
PB (E2615) - Recommended i	nded for less active		દે support	Privacy Flan	∩ 16" Width	Small Priva	ncv Flan)I			
PB Deep (E2615) - Recommer	nded for less active	e users seeking comfort &	È support		16" Width):			
PB Deep (E2615) - Recommer Lateral Elite Back Fi	ixed Laterals 5"Hx"	e users seeking comfort &	k support	Accessories: (○ 18" Width ○ 20" Width	Medium Pı Large Priva	ivacy Flap cy Flap	11			
PB Deep (E2615) - Recomment Lateral Elite Back Fi: Support Elite Back Sv	ixed Laterals 5"Hx' wing away Lateral	e users seeking comfort &		Accessories: (○ 18" Width ○ 20" Width	Medium Pı Large Priva	ivacy Flap				
PB Deep (E2615) - Recomment Lateral Elite Back Fit Support Elite Back Sv	ixed Laterals 5"Hx' wing away Lateral	e users seeking comfort & (7.5"L (E0956) Hardware (E0128))	Accessories: (○ 18" Width ○ 20" Width	Medium Pı Large Priva	ivacy Flap cy Flap				
PB Deep (E2615) - Recommer Lateral Elite Back Fir Support Elite Back Sv Swing away l	ixed Laterals 5"Hx wing away Lateral Lateral Pad (each)	7.5"L (E0956) Hardware (E0128) 4.25"Hx6.5"L (E0956) Mat	Accessories: (○ 18" Width ○ 20" Width	Medium Pı Large Priva	ivacy Flap cy Flap	11			
PB Deep (E2615) - Recomment Lateral Elite Back Fi: Support Elite Back Sv	ixed Laterals 5"Hx wing away Lateral Lateral Pad (each)	e users seeking comfort & (7.5"L (E0956) Hardware (E0128)	Mat	Accessories: (○ 18" Width ○ 20" Width	Medium Pı Large Priva	ivacy Flap cy Flap				
Description Processing	ixed Laterals 5"Hx' wing away Lateral Lateral Pad (each) © Elan Headre © Elan Standar © Elan 4 Point	ensers seeking comfort & (7.5"L (E0956) Hardware (E0128) 1 4.25"Hx6.5"L (E0956) est Hardware (E0128) and Pad with Standard (E04055)	Mat - Required with Select Cover (E0955)	Accessories: (18" Width 20" Width 22" Width	Medium Pr Large Priva Extra Large	ivacy Flap icy Flap e Privacy Flap				
Description Processing	ixed Laterals 5"Hx' wing away Lateral Lateral Pad (each) © Elan Headre © Elan Standar © Elan 4 Point	e users seeking comfort & (7.5"L (E0956) Hardware (E0128) 4.25"Hx6.5"L (E0956) est Hardware (E0128) ard Pad with Standard (Mat - Required with Select Cover (E0955)	Accessories: (18" Width 20" Width 22" Width	Medium Pi Large Priva Extra Large	ivacy Flap icy Flap e Privacy Flap				
Description Process Pr	ixed Laterals 5"Hx' wing away Lateral Lateral Pad (each) © Elan Headre © Elan Standar © Elan 4 Point	ensers seeking comfort & (7.5"L (E0956) Hardware (E0128) 1 4.25"Hx6.5"L (E0956) est Hardware (E0128) and Pad with Standard (E04055)	Mat - Required with Select Cover (E0955)	Accessories: (18" Width 20" Width 22" Width 22" Width	Medium Pr Large Priva Extra Large	ivacy Flap icy Flap e Privacy Flap				
Description Processing	ixed Laterals 5"Hx' wing away Lateral Lateral Pad (each) © Elan Headre © Elan Standar © Elan 4 Point	ensers seeking comfort & (7.5"L (E0956) Hardware (E0128) 1 4.25"Hx6.5"L (E0956) est Hardware (E0128) and Pad with Standard (E04055)	Mat - Required with Select Cover (E0955)	Accessories: (18" Width 20" Width 22" Width 22" Width	Medium Pr Large Priva Extra Large	ivacy Flap icy Flap e Privacy Flap				
Description Process Pr	ixed Laterals 5"Hx' wing away Lateral Lateral Pad (each) © Elan Headre © Elan Standar © Elan 4 Point	ensers seeking comfort & (7.5"L (E0956) Hardware (E0128) 1 4.25"Hx6.5"L (E0956) est Hardware (E0128) and Pad with Standard (E04055)	Mat - Required with Select Cover (E0955)	Accessories: (18" Width 20" Width 22" Width 22" Width	Medium Pr Large Priva Extra Large	ivacy Flap icy Flap e Privacy Flap				
PB Deep (E2615) - Recomment Lateral Elite Back Fire Elite Back Sw. Swing away be seen to be seen t	ixed Laterals 5"Hx' wing away Lateral Lateral Pad (each) Elan Headre Elan Standar Elan 4 Point Elan Headba	e users seeking comfort & (7.5"L (E0956) Hardware (E0128) 9 4.25"Hx6.5"L (E0956) est Hardware (E0128) ord Pad with Standard (t Pad (E0955)) and with 4 Point Pad (Mat - Required with Select Cover (E0955)	Accessories: (18" Width 20" Width 22" Width 22" Width	Medium Pr Large Priva Extra Large	ivacy Flap icy Flap e Privacy Flap				