

## Frequently Asked Questions / Answers

1. **What is the tip size on each curette?**
  - White FlexLoop is a 4mm wide tip.
  - Red AngleLoop is a 4mm wide tip.
  - Yellow CeraSpoon is a 4mm wide tip.
  - Green MicroLoop is a 3mm wide tip.
  - Purple VersaLoop is a 3mm wide tip.
  - Blue InfantScoop is a 2mm wide tip.
  - Orange ControlLoop is a 4mm wide tip.
2. **How do they come packaged?**

Fifty curettes per box.
3. **Are curettes individually wrapped?**

No, fifty curettes come in a single pouch for easy access.
4. **Are they sterile?**

Not necessary, not used in a sterile environment.
5. **What are the product numbers?**
  - # 9555 White FlexLoop
  - # 3444 Red AngleLoop
  - # 6333 Yellow CeraSpoon
  - # 1222 Green MicroLoop
  - # 4111 Purple VersaLoop
  - # 4888 Blue InfantScoop
  - # 2999 Orange ControlLoop
6. **Is there a CPT code for this procedure?**

We recommend that you submit this under CPT code 69210 for potential reimbursement. Alternatively, cerumen removal may increase the complexity of the visit. Curettage can also be used for foreign body removal. It can be submitted under CPT code 69200 for potential reimbursement.