

Positioning Sling



Positioning Sling for Full Body Support

- Provides full body support for lateral supine transfers, repositioning and turning in bed
- Multiple loop and strap options can accommodate various user repositioning and transfer tasks
- Intended to be left in place (assessment required)

Materials



Mesh



Poly




Spacer



Disposable

Applications



Supine transfer



Turning in bed



Higher up in bed

Part Numbers

Mesh

Part No.	Size	SWL
507810	36" x 70"	800 lb (363 kg)
507811	42" x 78"	800 lb (363 kg)
507814	56" x 78"	1,000 lb (454 kg)

Spacer with Elastic

Part No.	Size	SWL
9Y4901	42" x 78"	450 lb (204 kg)

Poly

Part No.	Size	SWL
507802	36" x 51"	800 lb (363 kg)
507807	36" x 70"	800 lb (363 kg)
507816	42" x 78"	800 lb (363 kg)
507818	56" x 78"	1,000 lb (454 kg)

Poly with Tab

Part No.	Size	SWL
507817	42" x 78"	800 lb (363 kg)

Disposable

Part No.	Size	SWL
547820	42" x 66"	1,100 lb (499 kg)
547810S	44" x 69"	660 lb (299 kg)
547810	52" x 78"	1,000 lb (454 kg)
547810XXL	63" x 87"	1,000 lb (454 kg)