

# BARISELECT®

Low Air Loss Mattress Replacement



## User-Friendly Control Unit

Rapid inflation, easy set-up and a motor so quiet, the patient will hardly know it's there!

**The SelectAir® System** Quietly circulates air through individual cells and allows air to escape from the mattress.

**Wicking Action:** Air escapes through tiny holes on the top of the mattress and circulates under a vapor-permeable cover to wick moisture away from patient.

## BARISELECT® Features

**Two extra inches** of mattress depth (8"), extended widths and specialized software accommodate extra weight capacity.

**Pulsation (alternating pressure):** Provides two cycle settings that continually readjust the system to ensure patient comfort.

**Manual Mode** allows you to input the patient's actual weight for optimal immersion.

**LCD Screen** with back-lit keys allowing system features to be visible in dark or dimly lit environments.

**Auto Comfort:** Takes the guess work out of programming. During set-up, mattress sets optimal pressure according to patient's weight distribution.

**Zoned, "true" low air loss used in acute care hospitals nationwide**

**Weight Limit:** 1,000 lbs. (453 kg)

**Mattress:** Fits most bariatric bed frames. Custom sizes available. Inquire for details.

**Cover:** Water resistant, vapor-permeable, flame resistant and biocompatible

**Warranty:** One Year Limited Warranty on control unit and soft goods

**Medicare Code:** E0277 Group 2 Support Surface

LIT-0005 | R100121

Manufactured in the USA by  
**MOXI**<sup>™</sup>  
ENTERPRISES  
Innovative health solutions