

ResQUp Product Specifications



Product Description	An innovative, lightweight, compact, and versatile mobility aid, therapeutic exercise tool, and fall recovery device designed to safely assist fallen persons up to a standing or seated position following a non-injurious fall.
Manufacturer	Safe Patient Solutions, LLC
Model Number	Model A
Color	Graystone (shown above)
Medical Device Classification	FDA Listed as a Class 1 Transfer Aid
Weight	15 lbs
Composition	Recyclable LDPE
Dimensions Folded	14" x 31" x 11" (H x W x D); Box Size: 15" x 31.5" x 11.5"
Dimensions Unfolded	14" x 31" x 28" (H x W x D)
Height of Each Level	4.5"
Weight Capacity	350 lbs
Fasteners	Zinc-plated steel
Non-Slip Pads	PORON® Non-absorbing Urethane
Hand Grips	Lateral Horizontal
USPTO Patent No.	8,393,684

Cleaning Recommendations

The ResQUp, like other mobile lifts, should be cleaned regularly. We recommend that the urethane pads and all composite surfaces on the ResQUp be cleansed using medical-grade antiseptic soap and water. Low-level disinfectants such as 70% - 90% alcohol, a 1:100 dilution of bleach, disinfectant wipes, and quaternary ammonium germicidal detergents may also be used. Do not use strongly acidic or basic solutions, phenol-based disinfectants, or organic solvents on any part of the ResQUp.