

## THERMAL VENT® OPEN KNEE WRAP STABILIZER

- Easy to apply and adjust
- Flexible stays help increase stability and support
- Opening behind knee helps increase comfort
- Heat therapy, support, and compression to help reduce pain
- Helps relieve aches, pains, strains, sprains, and weakened knees
- Bilateral design fits right or left knee



### THERMAL TECHNOLOGY

Engineered with MVT<sup>2</sup>:  
MicroVentilated Thermal Technology™

Thermal Vent Supports are a valuable aid to help rehabilitate injured tissue. These products provide warmth and heat therapy to help increase pliability and elasticity of muscles and tendons, making them more responsive to therapy and exercise.



### FUNCTIONAL & COMFORTABLE

The Swede-O® Thermal Vent® Open Knee Wrap Stabilizer provides comfortable support, gentle compression, and heat therapy to help relieve knee pain and promote healing. The flexible internal stays provide medial and lateral knee support for increased stability.

The wrap around design is easy and less painful to apply and adjust compared to pull-on style knee supports. The front opening reduces pressure across the patella which helps prevent aggravation to tenderness or bruising. Popliteal cut out helps increase comfort and reduces material bunching during knee flexion.

The adjustable closures provide a customized fit and allow for variations due to swelling. The material is soft and comfortable against the skin; bilateral design fits right or left knee.

### SPECIFICATIONS

#### Product Selection:

BRE-6453 Swede-O® Thermal Vent®  
Open Knee Wrap Stabilizer

**Product Weight:** .45 lb

**Individual Packaged Dimensions:**  
13" x 9" x 2"

**Materials:** Outer lining: Nylon, Mid Layer: Neoprene, Inner Lining: Nylon.

**Latex:** This product is not made with natural rubber latex.

MODEL #	SIZE	MEASUREMENT	
BRE-6453-GR-SML	SM	12 1/2-13 1/4"	32-34cm
BRE-6453-GR-MED	MED	13 1/4-14 1/2"	34-37cm
BRE-6453-GR-LRG	LRG	14 1/2-15 3/4"	37-40cm
BRE-6453-GR-1XL	XL	15 3/4-17"	40-43cm
BRE-6453-GR-2XL	2XL	17-18 1/4"	43-46cm
BRE-6453-GR-3XL	3XL	18 1/4-19 1/2"	46-50cm
BRE-6453-GR-4XL	4XL	19 1/2-21"	50-53cm
BRE-6453-GR-5XL	5XL	21-24"	53-61cm

Your Thermal Vent Support should be a snug and comfortable fit.  
Not too tight, not too loose. If in doubt, try next larger size.

Measure slightly bent  
beneath kneecap

