

## Introducing new ThermalSoft<sup>®</sup> Gel Hot and Cold Packs

ThermalSoft Gel packs are double-duty packs. Heat them in the microwave for heat therapy. For cold therapy put them in the freezer. ThermalSoft Gel Packs stay soft, even at freezer temperatures.

The large gel pieces mold to your body to deliver soothing heat or cold therapy. Two thin layers of vinyl make these packs sturdy yet conformable to the body.

A one year warranty against defects protects your investment.

And... your patients will love them. They can take them home and heat the packs in their microwave to receive the benefits of heat therapy. Or if you want them to use cold therapy, they can store the pack in the freezer until needed.

Contact your Mettler distributor today or call 1.800.854.9305 to find out how you can receive the soothing benefits of ThermalSoft Gel.

| P/N  | Description        |
|------|--------------------|
| 9831 | Standard Gel Pack  |
| 9832 | Half Size Gel Pack |
| 9833 | Cervical Gel Pack  |

**Size** 11" x 14" 7.5" x 11" 23" x 8" Contour

9834 Extra Large Gel Pack 11" x 21"

