

# Zing MPS (Multi-Position Stander) Zing Supine

## Size 1-Options

# Size 1

up to 44"  
(112cm)



### The MPS (multi-position stander) base includes:

- Four Swivel Locking Casters
- Foot Operated Gas Spring
- Mast Receiver with Prone Position Lockout

### The Supine base includes:

- Four Swivel Locking Casters
- Foot Operated Gas Spring
- Mast Receiver



PA5520 Zing MPS Base  
PA5522 Zing Supine Base

Comes with four swivel locking casters, foot operated gas spring, mast receiver. MPS base comes with prone position lockout.



PA5530 Mast

Mast is the base for positioning options including foot plates, knee supports, pelvic supports, upper body supports and head supports.



PA5532 Mast with Leg Abduction

Mast is the base for support and positioning components and allows for 0-30° of independent leg abduction. Multi-adjustable knee pads and multi-adjustable foot plates required for leg abduction.



PA5612 Dual Control

Allows both foot and hand control of the gas spring lift. Push handle PA5594 required.



PA5614 Gas Spring Lift Lockout

Locks gas spring in any position from lying down to fully upright, supine and prone positions.



PA5608 Directional Locking Caster

Great for transporting stander. Directionally locks one caster while leaving the other three free to swivel.



PA5536 Platform Foot Plate

Adjustable in height. Foot holders and foot straps can be added for more precise positioning of the feet.



PA5624, PA5626, PA5628  
Multi-Adjustable Foot Plates

Foot plates are adjustable forward/aft +/- .75" (2cm) +extra 1.5" (3.8cm), plantar/dorsi +/- 20°, toe-in/toe-out, swivel 30° left/right, and in height. Foot holder sizes: PA5624-6"Lx3"W (15x8cm), PA5626-7.75"Lx3.25"W (20x8cm), PA5628-9.75"Lx4"W (25x10cm).



PA5540, PA5542, PA5544  
Foot Holders for Platform Foot Plate

Use with the platform foot plate, gives added positioning of users feet. Foot holder sizes: PA5540-6"Lx3"W (15x8cm), PA5542-7.75"Lx3.25"W (20x8cm), PA5544-9.75"Lx4"W (25x10cm).



PA5546 Secure Foot Straps-8"L  
PT50082 Secure Foot Straps-10"L

Velcro® and D-ring adjustment holds feet securely in place at ankles and toes with multiple attachment slots for proper foot positioning. Straps measured in length over top of foot, with approximately 1.5" (4cm) of Velcro® overlap.



PA5548 Foot Straps-8"L  
PT50080 Foot Straps-10"L

Velcro® and D-ring adjustment holds feet in alignment with multiple attachment slots for proper foot positioning. Straps measured in length over top of foot, with approximately 1.5" (4cm) of Velcro® overlap.



PA5550 Calf/Knee Pad with Strap

4"x10" (10x25cm) pad attached to a 2" (5cm) Velcro® belt with a D-Ring closure, provides calf/knee stability for user. Fits circumference up to 24" (61cm). Gives support from the back of the calf/knee. Adjustable in height and fore/aft.



PA5556, PA5558, PA5729, PA5730  
Multi-Adjustable Knee Pads

Multi-adjustable knee pads are adjustable in swivel, height +/- .75" (2cm), PA5556 inside width 2.5" (6cm), PA5558 inside width 3.25" (8cm), PA5729 inside width 4.25" (11cm) and PA5730 inside width 5" (13cm). With a depth range of 5" (13cm).



PA5560 Multi-Adjustable Calf Pads 3.5"  
PA5678 Multi-Adjustable Calf Pads 4.5"

Great for adapting to knee contractures. Gives support from the back of the calf. Adjustable in height with knee pad and fore/aft 4" (10cm) independently. Must order multi-adjustable knee pads.



PA5562 Planar Pad-5"Hx7"W  
PA5564 Planar Pad-5"Hx9"W

Cushioned pad covered with Neoprene®. Pad sizes: PA5562-5"Hx7"W (13x18cm), PA5564-5"Hx9"W (13x23cm).



PA5566 Form to Fit Pad-5"Hx10"W  
PA5570 Form to Fit Pad-5"Hx12"W

Form to fit pad can be custom fit to the user. Pad sizes: PA5566-5"Hx10"W (13x25cm) PA5570-5"Hx12"W (13x30cm). Not recommended with leg abduction mast.



PA5574 Hip Supports

Keeps the hips in symmetrical alignment. Adjustable from 6"-10" (15-25cm) with fore/aft 4" (10cm). Pad size is 3.5"x3.5" (9x9cm). For planar pads only.



PA3011 Pelvic Positioning Belt

2" (5cm) wide belt has a Velcro® belt with a D-Ring closure, provides hip stability for user. Fits circumference up to 24" (61cm).



PA5610 Pelvic Padded Positioning Strap

5"x8" (13x20cm) pad attached to a 2" (5cm) belt, Velcro® belt with a D-Ring closure. Fits circumference up to 24" (61cm).



PA3028, PA3029 Hygienic Covers

Washable slip on cover available in five different patterns. For planar pads PA5562 and PA5564.



PA3030, PA3032 Hygienic Covers

Washable slip on cover available in five different patterns. For form to fit pads PA5566 and PA5570.



**PA5562 Planar Pad-5\"Hx7\"W**  
**PA5564 Planar Pad-5\"Hx9\"W**  
 Cushioned pad covered with Neoprene®.  
 Pad sizes: PA5562-5\"Hx7\"W (13x18cm),  
 PA5564-5\"Hx9\"W (13x23cm).



**PA5610**  
**Upper Body Padded Positioning Strap**  
 5\"x8\" (13x20cm) pad attached to a 2\"  
 (5cm) Velcro® belt with a D-Ring  
 closure. Fits circumference up to 24\"  
 (61cm).



**PA5600 Black Molded Swing-Away Tray-11\"Lx19\"W**  
 Swing-away tray is accessible in all  
 positions and provides anterior support  
 as user moves from supine to upright  
 and prone positions. Adjustable in  
 depth, height and tray angle. Tray size is  
 11\"Lx19\"W (28x48cm).



**PA5616 Planar Pad-7\"Hx7\"W**  
**PA5618 Planar Pad-7\"Hx9\"W**  
 Cushioned pad covered with Neoprene®.  
 Pad sizes: PA5616-7\"Hx7\"W (18x18cm),  
 PA5618-7\"Hx9\"W (18x23cm).



**PA3028, PA3029, PA3045, PA3046**  
**Hygienic Covers**  
 Washable slip on cover available in five  
 different patterns. For planar pads  
 PA5562, PA5564, PA5616, PA5618.



**PA5602 Clear Swing Away Tray-11\"Lx19\"W**  
 Swing-away tray is accessible in all  
 positions and provides anterior support  
 as user moves from supine to upright  
 and prone positions. Adjustable in  
 depth, height and tray angle. Tray size is  
 11\"Lx19\"W (28x48cm).



**PA5568 Form to Fit Pad-5\"Hx10\"W**  
**PA5620 Form to Fit Pad-7\"Hx10\"W**  
 Form to fit pad with built in lateral support that  
 can be custom fit to the user. Pad size for  
 PA5568 is 5\"Hx10\"W (13x25cm) and PA5620  
 is 7\"Hx10\"W (18x25cm).



**PA3031, PA3033, PA3047, PA3048**  
**Hygienic Covers**  
 Washable slip on cover available in five  
 different patterns. For form to fit pads  
 PA5568, PA5572, PA5620, PA5622.



**PA5604 Forearm Wings-Black Molded Tray**  
 Forearm wings give additional area to rest  
 arms when in standing position.



**PA5572 Form to Fit Pad-5\"Hx12\"W**  
**PA5622 Form to Fit Pad-7\"Hx12\"W**  
 Form to fit pad with built in lateral support that  
 can be custom fit to the user. Pad size for  
 PA5572 is 5\"Hx12\"W (13x30cm) and PA5622  
 is 7\"Hx12\"W (18x30cm).



**PA3039 Hygienic Covers for Lateral Supports**  
 Washable slip on cover available in five  
 different patterns. For lateral supports  
 PA5576.



**PA5606 Forearm Wings-Clear Tray**  
 Forearm wings give additional area to rest  
 arms when in standing position.



**PA5576 Lateral Supports**  
 Keeps the upper torso in symmetrical  
 alignment. Adjustable from 7-11\"  
 (18-28cm) with a 3.5\" (9cm) height range.  
 Pad size is 3.5\"Hx3.5\"W (9x9cm). For  
 planar pads only.



**PY5626 Head Support-5\"Hx8\"W**  
**PY5628 Head Support-6\"Hx10\"W**  
 Cushion size is 5\"Hx8\"W (13x20cm),  
 6\"Hx10\"W (15x25cm). Depth range is +1\"  
 to -2.75\" (+2 to -7cm) from backrest.  
 Maximum height from seat is 24\" (61cm).



**PA3038 Hygienic Padded Cover**  
 Washable slip on cover is padded and  
 available in five different patterns.



**PA5578 Lateral Supports with Elbow Stop**  
 Works as a lateral support and elbow rest.  
 Great for the supine position. For planar  
 pads only.



**PA5586 Head Support-Form to Fit**  
 Cushion size 5\"Hx10\"W (13x25cm). Depth  
 range is -2.5\" (6cm) to +1.5\" (4cm) from  
 backrest. Maximum height from top of  
 upper body pad is 8\" (20cm).



**PNG50042 Hand Grips**  
 Provides additional arm and hand support  
 and alignment. Easily movable for precise  
 arm positioning. Attaches with suction  
 cups. Clear tray required.



**PA5580**  
**Lateral Supports w/Elbow Stop/Arm Rest**  
 Works as a lateral support, elbow rest  
 and armrest. Great for the supine, upright  
 and prone position. For planar pads only.  
 Height and width adjustable.



**PA3035, PA3036 Hygienic Covers**  
 Washable slip on cover available in five  
 different patterns. For PY5626 and  
 PY5628.



**PA5594 Push Handle**  
 Makes pushing the stander easier for the  
 caregiver and also aids in lifting and  
 lowering the child. Required for  
 dual control.



**PA5630 Elbow Stop with Arm Rest**  
 Great for loading and bringing to upright.  
 Works with both planar and form to fit  
 upper body supports. For planar and form  
 to fit pads. Height and width adjustable.



**PA3037 Hygienic Cover**  
 Washable slip on cover available in five  
 different patterns. For PA5586.



**P80809 Tool Pouch**  
 Holds the assembly tools provided with  
 the stander. Attaches to unit with Velcro®  
 straps.



**PT50086 X-Style Chest Vest-9.5\"Lx9\"W**  
**PT50088 X-Style Chest Vest-11\"Lx9.5\"W**  
 Neoprene® with camlock adjustment  
 provides anterior trunk support for child.



**PA5596 Black Molded Tray-11\"Lx19\"W**  
 Tray is accessible in all positions and  
 provides anterior support as user moves  
 from supine to upright and prone  
 positions. Adjustable in depth and height.  
 Tray size is 11\"Lx19\"W (28x48cm).



**P80246 Angle Locator**  
 Magnetic angle locator can be utilized for  
 standing progress and leg abduction.

# Zing MPS (Multi-Position Stander) Zing Supine

## Size 2-Options

# Size 2

40"-60"  
(102-152cm)



### The MPS (multi-position stander) base includes:

- Four Swivel Locking Casters
- Foot Operated Gas Spring or Pow'r Up Lift
- Mast Receiver with Prone Position Lockout

### The Supine base includes:

- Four Swivel Locking Casters
- Foot Operated Gas Spring or Pow'r Up Lift
- Mast Receiver



**PB5500 Zing MPS Base Gas Spring Lift**  
**PB5502 Zing Supine Base Gas Spring Lift**  
Comes with four swivel locking casters, foot operated gas spring, mast receiver. MPS base comes with prone position lockout.



**PB5504 Zing MPS Base Pow'r Up Lift**  
**PB5506 Zing Supine Base Pow'r Up Lift**  
Comes with four swivel locking casters, Pow'r Up Lift, mast receiver. MPS base comes with prone position lockout.



**PB5512 Mast**  
Mast is the base for positioning options including foot plates, knee supports, pelvic supports, upper body supports and head supports.



**PB5514 Mast with Leg Abduction**  
Mast is the base for support and positioning components and allows for 0-30° of independent leg abduction. Multi-adjustable knee pads and multi-adjustable foot plates required for leg abduction.



**PB5516 Dual Control**  
Allows both foot and hand control of the gas spring lift. Push handle PB5560 required.



**PA5614 Gas Spring Lift Lockout**  
Locks gas spring in any position from lying down to fully upright, supine and prone positions.



**PB5510 Directional Locking Caster**  
Great for transporting stander. Directionally locks one caster while leaving the other three free to swivel.



**PB5518 Platform Foot Plate**  
Adjustable in height. Foot holders and foot straps can be added for more precise positioning of the feet. Platform 16.75"Wx11.5"L (29-43cm).



**PB5520, PB5522, PB5524 Multi-Adjustable Foot Plates**  
Foot plates are adjustable forward/aft +/- .75" (2cm) +extra 1.5" (3.8cm), plantar/dorsi +/- 20°, toe-in/toe-out, swivel 30° left/right, and in height. Foot holder sizes: PB5520 7.75"Lx3.25"W (20x8cm), PB5522 9.75"Lx4"W (25x10cm), PB5524 11.75"Lx5"W (30x13cm).



**PA5542, PA5544, PB5526 Foot Holders for Platform Foot Plate**  
Use with the platform foot plate, provides added positioning of users feet. Foot holder sizes: PA5542 7.75"Lx3.25"W (20x8cm), PA5544 9.75"Lx4"W (25x10cm), PB5526 11.75"Lx5"W (30x13cm).



**PT50082 Secure Foot Straps-10"L**  
**PY5638 Secure Foot Straps-15"L**  
**PNG50506 Secure Foot Straps-21"L**  
Velcro® and D-ring adjustment holds feet securely in place at ankles and toes with multiple attachment slots for proper foot positioning. Straps measured in length over top of foot, with approximately 1.5" (4cm) of Velcro® overlap.



**PT50080 Foot Straps-10"L**  
**PY5636 Foot Straps-15"L**  
**PNG50505 Foot Straps-21"L**  
Velcro® and D-ring adjustment holds feet in alignment with multiple attachment slots for proper foot positioning. Straps measured in length over top of foot, with approximately 1.5" (4cm) of Velcro® overlap.



**PB5528 Calf/Knee Pad with Strap**  
17.5"x5.5" (44x14cm) pad attached to a 2" (5cm) Velcro® belt with a D-Ring closure, provides calf/knee stability for user. Fits circumference up to 28"-38" (71-97cm). Gives support from the back of the calf/knee. Adjustable in height and fore/aft.



**PB5530, PB5532, PB5534 Multi-Adjustable Knee Pads**  
Multi-adjustable knee pads are adjustable in swivel, height +/- .75" (2cm), PA5530 inside width 3.25" (8cm), PA5532 inside width 4.25" (11cm) and PB5534 inside width 5" (13cm). With a depth range of 7" (18cm).



**PA5678 Multi-Adjustable Calf Pads 4.5"**  
**PB5536 Multi-Adjustable Calf Pads 5.5"**  
Great for adapting to knee contractures. Gives support from the back of the calf. Adjustable in height with knee pad and fore/aft 4" (10cm) independently. Must order multi-adjustable knee pads.



**PB3000 Planar Pad 10"Hx11"W**  
**PB3002 Planar Pad 10"Hx13"W**  
Cushioned pad covered with Neoprene®. Pad sizes: PB3000-10"Hx11"W (25x28cm) PB3002-10"Hx13"W (25x33cm)



**PB3004 Form to Fit Pad 10"Hx11"W**  
**PB3006 Form to Fit Pad 10"Hx13"W**  
Form to fit pad can be custom fit to the user. Pad sizes: PB3004 10"Hx11"W (25x28cm) PB3006 10"Hx13"W (25x33cm). Not recommended with leg abduction mast.



**PB5538, PB5540 Hip Supports**  
Keeps the hips in symmetrical alignment. Adjustable from PB5538 7"-13" (18-33cm) and PB5540 12"-18" (30-46cm) with fore/aft 4" (10cm). Pad size is 4.5"x4" (11x10cm). For planar pads only.



**PT30066 Pelvic Positioning Strap**  
2" (5cm) wide belt has a Velcro® belt with a D-Ring closure, provides hip stability for user. Fits circumference up to 35" (89cm).



**PB5542 Pelvic Padded Positioning Strap**  
6"x10" (15x25cm) pad attached to a 2" (5cm) Velcro® belt with a D-Ring closure. Fits circumference up to 35" (89cm).



**PB3025, PB3026 Hygienic Cover**  
Washable slip on cover available in five different patterns. For planar pads PB3000 and PB3002.

**PB3027, PB3028 Hygienic Cover**

Washable slip on cover available in five different patterns. For form to fit pads PB3004 and PB3006.

**PB5542 Upper Body Padded Positioning Strap**

6"x10" (15x25cm) pad attached to a 2" (5cm) Velcro® belt with a D-Ring closure. Fits circumference up to 35" (89cm).

**PB5558 Clear Swing Away Tray-21"Lx24"W**

Swing-away tray is accessible in all positions and provides anterior support as user moves from supine to upright and prone positions. Adjustable in depth, height and tray angle. Tray size is 21"Lx24"W (53x61cm).

**PB3001 Planar Pad-8"Hx11"W  
PB3003 Planar Pad-8"Hx13"W**

Cushioned pad covered with Neoprene®. Pad size for PB3001 is 8"Hx11"W (20x28cm) and PB3003 is 8"Hx13"W (20x33cm).

**PB3029, PB3030 Hygienic Cover**

Washable slip on cover available in five different patterns. For planar pads PB3301 and PB3003.

**PT50156 Elbow Pad-9.25" Cutout**

Additional padded support surface on the black molded or clear tray. 9.25" (24cm) cutout width. Overall width 21.5" (55cm), pad extends back 3.75" (10cm) posterior from current table edge for additional upper body/extremity support.

**PB3005 Form to Fit Pad-8"Hx11"W  
PB3007 Form to Fit Pad-8"Hx13"W**

Form to fit pad with built in lateral support that can be custom fit to the user. Pad size for PB3005 is 8"Hx11"W (20x28cm) and PB3007 is 8"Hx13"W (20x33cm).

**PB3031, PB3032 Hygienic Covers**

Washable slip on cover available in five different patterns. For form to fit pads PB3005 and PB3007.

**PT50160 Elbow Pad-11" Cutout**

Additional padded support surface on the black molded or clear tray. 11" (28cm) cutout width. Overall width 21.5" (55cm), pad extends back 3.75" (10cm) posterior from current table edge for additional upper body/extremity support.

**PB5544, PB5546 Lateral Supports**

Keeps the upper torso in symmetrical alignment. Adjustable from PB5544 7"-11"W with PB3001 and 9"-13"W with PB3003 and PB5546 10"-14"W with PB3001 and 12"-16"W with PB3003 with a 3.5" (9cm) height range. Pad size is 4"Hx6"W (10x15cm). For planar pads only.

**PB3049 Hygienic Cover for PB5544/PB5546**

Washable slip on cover available in five different patterns. For PB5544 and PB5546.

**PY5582 Elbow Pad-13" Cutout**

Additional padded support surface on the black molded or clear tray. 13" (33cm) cutout width. Overall width 21.5" (55cm), pad extends back 3.75" (10cm) posterior from current table edge for additional upper body/extremity support.

**PB5548 Lateral Supports with Elbow Stop**

Works as a lateral support and elbow rest. Great for the supine position. Range: 10"-17" (25-43cm) with planar pad B3001 and 12"-19" (30-48cm) with planar pad B3003. For planar pads only.

**PY5626 Head Support 5"Hx8"W  
PY5628 Head Support 6"Hx10"W  
PY5630 Head Support 7"Hx14"W**

Cushion sizes are: PY5626 5"Hx8"W (13x20cm), PY5628 6"Hx10"W (15x25cm), PY5630 7"Hx14"W (18x36cm). Depth range is +1" to -2.75" (+2 to -7cm) from backrest. Maximum height from top of upper body pad is 10" (25cm).

**PNG30295 Hygienic Padded Cover**

Washable slip on cover is padded and available in five different patterns.

**PB5550 Lateral Supports w/Elbow Stop/Arm Rest**

Works as a lateral support, elbow rest and armrest. Great for the supine, upright and prone position. For planar pads only. Height and width adjustable.

**PA5586 Head Support-Form to Fit**

Cushion size 5"Hx10"W (13x25cm). Depth range is -2.5" (6cm) to +1.5" (4cm) from backrest. Maximum height from top of upper body pad is 10" (25cm).

**PNG50042 Hand Grips**

Provides additional arm and hand support and alignment. Easily movable for precise arm positioning. Attaches with suction cups. Clear tray required.

**PB5552 Elbow Stop with Arm Rest**

Great for loading and bringing to upright. Works with both planar and form to fit upper body supports. For planar pads only. Height and width adjustable.

**PA3035, PA3036, PB3033, PA3037 Hygienic Cover**

Washable slip on cover available in five different patterns. For PY5626, PY5628 PY5630 and PA5586.

**PB5560 Push Handle**

Makes pushing the stander easier for the caregiver and also aids in lifting and lowering the child. Required for dual control.

**PNG30170 X-Style Chest Vest-9.5"Lx9"W  
PNG30056 X-Style Chest Vest-11"Lx9.5"W**

Neoprene® with camlock adjustment provides anterior trunk support for child.

**PB5554 Black Molded Tray-21"Lx24"W**

Tray is accessible in all positions and provides anterior support as user moves from supine to upright and prone positions. Adjustable in depth and height. Tray size is 21"Lx24"W (53x61cm).

**P80809 Tool Pouch**

Holds the assembly tools provided with the stander. Attaches to unit with Velcro® straps.

**PNG30171 Y-Style Chest Vest-17.5"Lx11.5"W  
PNG30021 Y-Style Chest Vest-20.5"Lx12.5"W**

Neoprene® with camlock adjustment provides anterior trunk support for child.

**PB5556 Black Molded Swing-Away Tray-21"Lx24"W**

Swing-away tray is accessible in all positions and provides anterior support as user moves from supine to upright and prone positions. Adjustable in depth, height and tray angle. Tray size is 21"Lx24"W (53x61cm).

**P80246 Angle Locator**

Magnetic angle locator can be utilized for standing progress and leg abduction.